

END-TO-END
COMPARISON CHART



Weight Capacity
Seat Width
Seat Depth
Seat-to-Floor
Chair Width
Chair Depth
Chair Height
Operation
Fabric Options

Alperton	Crosby	Ealing	Eton	Hoxton	Ludlow	Marcos	Mercer	Seagrove	Studio	Kennington
Dual motor lift lay-flat recliner with headrest and lumbar	Three motor lay-flat recliner with headrest and lumbar	Dual motor lift lay-flat recliner with headrest and lumbar	‘Box in box’ dual motor lift chair with headrest	Single motor lift chair	Dual motor lift chair with headrest and lumbar	Dual motor lift lay-flat recliner with headrest, leg extension and lumbar	Dual motor lift chair with headrest and lumbar	Dual motor lift lay-flat recliner with headrest and lumbar	Dual motor lift chair with headrest and lumbar	Duo seat, dual motor lift lay-flat recliner with headrest and lumbar
158.7 kg	136 kg	158.7 kg	180 kg	113.3 kg	158.7 kg	158.7 kg	158.7 kg	158.7 kg	158.7 kg	158.7 kg
43.1 cm	53.3 cm	53.3 cm	53.3 cm	49.5 cm	53.3 cm	57.1 cm	50.8 cm	55.8 cm	57.1 cm	57.1 cm
46.9 cm	53.3 cm	48.2 cm	52 cm	50.8 cm	53.3 cm	55.8 cm	54.6 cm	53.3 cm	54.6 cm	54.6 cm
45.7 cm	48.2 cm	49.5 cm	48.2 cm	49.5 cm	46.9 cm	49.5 cm	46.9 cm	49.5 cm	46.9 cm	49.5 cm
74.9 cm	85 cm	78.7 cm	87.6 cm	82.5 cm	83.8 cm	88.9 cm	86.4 cm	99 cm	96.5 cm	1.93 m
93.9 cm	105.5 cm	101.6 cm	95.3 cm	99 cm	104.2 cm	105.4 cm	104 cm	107.9 cm	106.7 cm	106.6 cm
99 cm	106.7 cm	106.6 cm	106.7 cm	101.5 cm	105.5 cm	107.9 cm	106.7 cm	106.6 cm	106.7 cm	106.6 cm
USB Hand Wand	Side Button/ USB Hand Wand	USB Hand Wand	USB Hand Wand	USB Hand Wand	USB Hand Wand	USB Hand Wand	USB Hand Wand	USB Hand Wand	USB Hand Wand	USB Hand Wand
Merino (Fresh) (Latte)	Stonewash (Gun Metal)	Radiance (Chestnut) Stark (Onyx)	Stonewash (Amber) (Dune)	Stonewash (Glove) (Mossy)	Stonewash (Dove) (Fudge)	Mosaic (Fudge) (Charcoal)	Badlands (Saddle) (Eclipse)	Picasso (Walnut) (Midnight)	Santa Barbara (Latte) (Bark) (Fog)	Stonewash (Gun Metal)