

## TurnTable

Used for transferring a user between a wheelchair and chair, toilet or bed. The TurnTable is placed on the floor under the user's feet allowing them to move. It is an outstanding, easy-to-use aid for skilled caregivers for transferring users.



### Specifications

Product specifications	
SWL	135kg
Diameter	30/38mm
Warranty	12 months

### Features

- The surface provides a good grip against both feet and floor
- The two plates rotate smoothly over each other
- A practical handle facilitates handling
- Convenient size and low weight
- Easy to clean