

# xtaga Impala

The ultimate in portability, power and performance. The xtaga Impala is packed with features and perfect for those seeking a compact scooter without having to compromise on luxury or performance.



## Colours Available



Puya Teal



Cobalt Blue

## Features



**Seat**  
Mid back captain seat with head rest



**Weight Capacity**  
145kg



**Wheels**  
4 Wheel



**Lighting**  
LED headlight, front basket, reflectors, brake lights

# Full Features

- Quick disassembly for transporting
- Height and width adjustable comfort padded armrests
- Pneumatic tyres for excellent outdoor performance
- Easy drive adjustable tiller with ergonomic wraparound handles
- Contoured, mid-back captain seat with padded headrest
- Large weight capacity
- Front & rear suspension
- Front basket
- Seat pocket for storage
- LED lower headlight & reflectors
- Horn
- Sleek, sporty style

# Specifications

| SPECIFICATIONS                | XTAGA IMPALA                            |
|-------------------------------|---|
| Maximum Weight Capacity       | 145kg                                   |
| Maximum Speed <sup>1</sup>    | 8km/h                                   |
| Maximum Range <sup>1</sup>    | Up to 29km                              |
| Length <sup>2</sup>           | 112cm                                   |
| Width <sup>2</sup>            | 56cm                                    |
| Scooter Weight <sup>2,3</sup> | 59.5kg / 82.5kg with Battery            |
| Ground Clearance <sup>2</sup> | 9cm                                     |
| Seat to Ground                | 47cm (Adjustable)                       |
| Turning Radius <sup>2</sup>   | 134.5cm                                 |
| Battery Requirements          | 50Ah                                    |
| Charger                       | Input AC 100-240V, Output DC24V, 3.5Amp |
| Seating                       | Captain Seat with Headrest              |
| Motor                         | 420W                                    |

# Compliance

| COMPLIANCE                                      |  |
|---|--|
| Australian Register of Therapeutic Goods (ARTG) | xtaga Impala is entered in the ARTG 406770 |
| ISO   | ISO9001, ISO13485                          |

<sup>1</sup>Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition and tyre condition. These specifications can be subject to a variance of (+/- 10%).

<sup>2</sup>Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of (+ or - 3%).

<sup>3</sup>Battery weight may vary based on manufacturer.

**IMPORTANT :** Always be sure to read your owner's manual thoroughly before operating your machine.

**WARNING :** Do not exceed 10km/ph on public roads or paths.

## Silverwater

Unit 8, 110-120 Silverwater Rd  
Silverwater NSW 2128  
Ph: 9649 2111

## Newcastle

7 Poynton Pl  
Thornton, NSW 2322  
Ph: 4047 3608