


Medical Scales



Active Mobility



Baby Scales



Hoist Scales



Wheelchair Scales



Bariatric Scales

Portable Baby Scale 3552011

Features

- Lightweight (3.1kg)
- Portable, ideal for health visitors and midwives
- Optional rucksack or carry case for easy transportation
- Large, clear display
- Suitable for neonatal baby units within hospitals (2g accuracy)
- Extremely long battery life (>600hrs.)



Optional Trolley, Carry Case and/or Rucksack

Specifications

Capacity	15kg
Graduations	2g<6kg>5g
Dimensions	370mm x 590mm
Weight of scale	3.1kg



Portable Baby Scale 3552004



Features

- Class III approved
- Compact, stable and robust
- Lightweight
- Secure hygienic slide on/off weighing pan
- Ideal for health visitors and midwives
- Ideal for weighing baby or stand-on toddler
- MDD approved

Specifications

Capacity	20kg
Graduations	5g<10kg>10g
Base dimensions	320mm x 310mm x 60mm
Tray dimensions	585mm x 280mm x 120mm
Weight of scale	4kg



Optional Trolley, Carry Case and Rucksack



Portable Baby Scale 3552006



Features

- Class III approved
- Compact, stable and robust
- Lightweight
- Secure hygienic slide on/off weighing pan
- Ideal for health visitors and midwives
- Ideal for weighing school children
- Optional rucksack or carry case
- MDD approved

Specifications

Capacity	50kg
Graduations	10g<20kg>20g
Base Dimensions	320mm x 310mm x 60mm
Tray Dimensions	585mm x 280mm x 120mm
Weight of scale	3.3kg



Optional Trolley, Carry Case and Rucksack



Lightweight Portable Baby Scale 3552003



Features

- Class III approved
- Ultra lightweight
- Secure hygienic slide on/off weighing pan
- Ideal for health visitors and midwives
- Optional rucksack or carry case

Specifications

Capacity	20kg
Graduations	10g<10kg>20g
Base dimensions	560mm x 290mm x 120mm
Tray dimensions	560mm x 290mm
Weight of scale	2.4kg



Optional Trolley, Carry Case or Rucksack



Hand Held Baby Scale 3552008

Features

- Ultra lightweight
- Portable
- Ideal for health visitors and midwives
- Reusable sling included
- Optional carry case

Specifications

Capacity	10kg
Graduations	10g
Dimensions	79mm x 150mm x 193mm
Weight of scale	0.45kg



3 in 1 - Baby, Toddler, Adult Scale 3552002

Features

- Ideal for mother and baby clinics and GP Practices
- Removable baby bowl
- Seat for toddlers
- Large base for adults
- Heavy duty construction
- Automatic selection of correct weighing range
- Optional printer

Specifications

Capacity	150kg
Graduations	Baby/Toddler: 60kg x 20g Auto Select Multi-range Adult: 150kg x 50g Auto Select
Platform Dimensions	700mm (w) x 600mm (d) x 45mm (h)
Adapter (supplied)	3pin tip +ve 12V
Weight of scale	37kg



Digital Portable Scale 3552016

Features

- Ultra lightweight
- Extremely accurate
- Reversed weight reading option for discreet reading
- Ideal for home visits
- Optional carry case for easy transportation
- MDD approved
- **OPTIONAL BLUETOOTH (3552017)**

Specifications

Capacity	220kg
Graduations	50g<150kg<100g
Dimensions	3pin tip +ve 12V
Weight of scale	3kg



Optional Carry Case



Portable Scale with Remote Display 3552018



Features

- Remote display
- Ultra lightweight
- Extremely accurate
- Reversed weight reading option for discreet reading
- Ideal for home visits
- Optional carry case for easy transportation
- MDD approved

Specifications

Capacity	220kg
Graduations	50g<150kg>100g
Dimensions	320mm x 310mm x 60mm
Weight of scale	3kg



Optional Carry Case



Digital Portable Scale 3552020



Features

- Lightweight
- Large LCD display with reversible read function
- Non-slip surface
- Ideal for home visits
- Optional carry case for easy transportation
- MDD approved
- **OPTIONAL BLUETOOTH (3552021)**

Specifications

Capacity	220kg
Graduations	200g
Dimensions	320mm x 310mm x 60mm
Weight of scale	3kg




Optional Carry Case
6 | www.activemobility.com.au



Portable Floor Scale with Remote indicator 3552023

Features

- Lightweight and robust
- Suitable for hospitals, GP surgeries and health visitors
- Optional carry case for easy transportation
- MDD approved
- **HINGED WEIGHT INDICATOR OPTION (3552024)**



Specifications

Capacity	250kg
Graduations	100g
Dimensions	310mm (w) x 90mm(h) x 310mm (d)
Weight of scale	6kg



Optional Carry Case | Hinged weight indicator model (3552024)



Portable Floor Scale 3552027

Features

- Extra wide, low profile platform suitable for bariatric and geriatric patients
- Discreet weighing
- Light platform
- High capacity
- Optional carry case
- Ideal for GP surgeries or hospitals
- **500kg MODEL OPTION (3552028)**



Specifications

Capacity	300kg (3552027) / 500kg (3552028)
Graduations	50g<150kg>100g (300kg) 100g<300kg>200g (500kg)
Dimensions	600mm (w) x 50mm (h) x 390mm (d)
Weight of scale	8kg




Optional Carry Case



Column Scale with Integrated Height Measure 3552040

Features

- Unique design
- Heavy duty construction
- Large platform to suit obese patients
- Integrated height measure
- High capacity
- Suitable for outpatient and bariatric clinics
- Optional printer



Specifications

Capacity	300kg
Graduations	50g<150kg>100g
Dimensions	440mm (w) x 1240mm (h) x 460mm (d)
Weight of scale	15kg



Professional Column Scale 3552041



Features

- High capacity
- Excellent accuracy
- Ideal for outpatient clinics
- Optional digital or manual height measures
- Optional printer
- MDD approved

Specifications

Capacity	250kg
Graduations	100g
Base Dimensions	310mm (w) x 310mm (d) x 90mm (h)
Weight of scale	8kg



Optional Thermal Printer



HOLD

TARE



MDD

BMI

Low Cost Column Scale 3552045



Features

- Low cost, high capacity scale
- Light and easy to move
- Suitable for GP surgeries and outpatient clinics
- Optional height measure

Specifications

Capacity	250kg
Graduations	200g<125kg>500g
Dimensions	310mm (w) x 90mm (h) x 310mm (d)
Weight of scale	8kg



HOLD

TARE



BMI



Optional height measure

High Capacity Chair 3552000

Features

- Extra wide super-sized seat
- Robust construction
- Four wheels (two with brake)
- Hinged footrests for patients comfort
- Optional seat belt or harness
- Optional printer
- MDD approved
- 300kg CHAIR SCALE AVAILABLE (3552001)



Specifications

Capacity	250kg / 300kg
Graduations	50g<150kg>100g
Beam Dimensions	530mm (w) x 460mm (d)
Weight of scale	27kg



Option:
No footrests



Option:
Fold away to the side footrest

Professional Chair Scale 3552047



Features

- Light and easy to move
- Four wheels (two with brake)
- Hinged armrests for easier patient transfer
- Hinged footrests for patient comfort
- Optional seat belt or harness
- Optional printer



Specifications

Capacity	250kg
Graduations	50g<150kg>100g
Seat Dimensions	400mm (w) x 380mm (d)
Weight of scale	18kg



Low Cost Chair Scale 3552048



Features

- Light and easy to move
- Four wheels (two with brake)
- Hinged armrests for easier patient transfer
- Hinged footrests for patient comfort
- Optional seat belt or harness

Specifications

Capacity	250kg
Graduations	200g<125kg>500g
Seat Dimensions	400mm (w) x 380mm (d)
Weight of scale	18kg



Professional Wheelchair Scale 3552053

Features

- One-piece and multifunctional
- Large weighing area
- Fitted castors for easy transportation
- Easy to store
- 2 ramps
- Optional column for indicator
- Optional handrails
- Tare to deduct weight of wheelchair
- Optional printer



Specifications

Capacity	300kg
Graduations	100g
Beam Dimensions	1150mm (w) x 800mm (d) x 66mm (d)
Weight of scale	28.6kg



Optional column



Optional Handrails

Professional Wheelchair Scale 3552049



Features

- Low profile platform
- Tare to deduct weight of the wheelchair
- Folded for easy transport and storage
- Easy access with roll on / roll off side
- Optional printer



Specifications

Capacity	300kg
Graduations	100g
Beam Dimensions	700mm (w) x 800mm (h) x 80mm (d)

CLASS
III
APPROVED

HOLD

TARE



BMI

Portable Wheelchair Weigh Beams 3552050



Features

- Two portable weigh beams
- Lightweight and portable for use in the community
- Tare to deduct weight of the wheelchair
- Easy storage
- Optional carry case
- Optional printer



Specifications

Capacity	300kg
Graduations	100g
Beam Dimensions	1200mm (l) x 150mm (w) x 95mm (h)
Weight of scale	8.5kg



Optional Carry Case

CLASS
III
APPROVED

HOLD

TARE




BMI

MDD

Floor Mounted Bed Scale with Ramps 3552100

Features

- Made to measure weighing area
- Capacity and graduation to order
- Anti slip durbar finish
- 2 long ramps giving shallow incline
- Tare to deduct weight of bed
- Wall mounted or column mounted indicator
- Optional printer



Specifications


Capacity	600kg to 1500kg
Graduations	200g to 500g
Base Dimensions	1200mm x 1000mm 1200mm x 1200mm 1500mm x 1500mm (other dimensions available on request)



Stainless Steel Bed Weighing Scale with Inbuilt Ramps 3552101

Features

- Drive thru platform scale with built in foldable ramps
- Very low profile platform only 45mm high
- Short ramps with shallow incline
- Fully stainless steel construction
- Tare to deduct weight of bed
- Optional printer
- Wall mounted or column mounted indicator



Specifications


Capacity	600kg / 1500kg
Graduations	200g / 500g
Dimensions	1200mm x 1200mm 1200mm x 1500mm 1500mm x 1500mm



Pit Mounted Bed Scale 3552102

Features

- Pit mounted platform (top of scale is level to floor)
- Capacity and graduation to order
- No need for ramps
- Anti slip durbar finish
- Easy loadcell access
- Tare to deduct weight of bed
- Optional printer
- Wall mounted or column mounted indicator



Specifications

Capacity	600kg to 1500kg
Graduations	200g to 500g
Base dimensions	1200mm x 1200mm 1500mm x 1500mm Other dimensions available on request
Power supply	3pin tip +ve 12V



Bed Scale 3552098

Features

- Easy and accurate weighing of bedridden patients
- Portable bed weighing beams
- Tare to deduct weight of the bed
- Easy to handle and store
- Optional printer



Specifications

Capacity	600kg
Graduations	200g
Dimensions	880mm (l) x 120mm (w) x 15mm (h)
Weight of beam (each)	11.5kg



Carry Handle



High Capacity Handrail Scale 3552036



Features

- Ultra stable and robust
- Flat packed for easy transport
- Wheels for easy movement
- Large, non-slip platform
- Suitable for bariatric and geriatric patients
- Optional printer



Specifications

Capacity	300kg
Graduations	50g<150kg>100g
Platform dimensions	550mm (l) x 550mm (w) x 1160mm (h)
Weight of scale	17.7kg

CLASS
III
APPROVED

HOLD

TARE



BMI

MDD

Professional Personal and Wheelchair Scale 3552052



Features

- Large platform with built-in ramp for wheelchair access
- Hinged seat for less stable patients
- Handrail to support bariatric and geriatric patients
- Wheels for easy movement
- Contoured ramp for easy access from two sides
- Optional printer



Specifications

Capacity	300kg or 500kg
Graduations	300kg: 100g 500kg: 100g<300kg>200g
Platform dimensions	800mm x 800mm
Weight of scale	46kg

CLASS
III
APPROVED

HOLD

TARE



BMI

High Capacity Handrail Scale 3552035



Features

- Curved, strong handrails for supporting patient
- Low-level platform for easy access
- Wheels for easy movement
- Large platform suitable for bariatric and geriatric patients
- Indicator swivels for optimum viewing angle
- Optional digital or manual height measure
- Optional printer

Specifications

Capacity	300kg
Graduations	50g<150kg>100g
Platform dimensions	700mm (w) x 600mm (d) x 50mm (h)
Weight of scale	33kg



CLASS
III
APPROVED

HOLD

TARE



BMI

Hoist Weighing Attachment 3552056



Features

- Lightweight and portable for use in care homes, community or hospital
- Suits most hoists 'coathanger' sling support
- Safety clips to secure sling straps
- Optional carry case
- MDD approved

Specifications

Capacity	200kg
Graduations	100g
Dimensions	630mm (l) x 275mm (h) x 160mm (d)
Weight of scale	6kg



Carry Case



Hoist Weighing Attachment 3552089



Features

TO SUIT LARGER HOISTS

- Lightweight and portable for use in home, community or hospital
- Suits most hoists 'coathanger' sling support
- Safety clips to secure sling straps
- Larger design to suit larger hoists
- Optional carry case
- MDD approved

Specifications

Capacity	200kg
Graduations	100g
Dimensions	730mm (l) x 275mm (h) x 160mm (d)
Weight of scale	6.7kg



Carry Case



Hoist Weighing Module 2211065

Features

- Universal weighing module with tare function and kg/lb conversion
- Suitable for yoke or pivot attachments
- Available on all Allegro brand hoists and some non-Allegro hoists (as approved by AMS)
- Useful where weighing of patients is difficult due to the level of support required



Specifications

Capacity	300kg
Graduations	100g
Dimensions	95mm (w) x 170mm (h) x 80mm (d)
Weight of scale	1kg

Hoist Weigh Module / Integrated Yoke 2211037

Features

- Reduces loss of lifting height compared with other hoist scales
- Use with general purpose (loop attachment) slings
- Non-absorbent polyurethane padding



Specifications

Capacity	300kg
Graduations	100g
Dimensions	570mm (w) x 200mm (h) x 105mm (d)
Weight of scale	1kg

Allegro Hoist Weigh Module / Integrated Pivot Frame 2211039



Features

- Reduces loss of lifting height compared with other hoist scales
- Use with pivot (clip attachment) slings
- Non-absorbent polyurethane padding



Specifications

Capacity	300kg
Graduations	100g
Dimensions	550mm (w) x variable height and length
Weight of scale	8kg

Allegro Hoist Weigh Module / Integrated Power Pivot Frame 2211041



Features

- Improved manual handling
- Reduces loss of lifting height compared with other hoist scales
- Use with pivot (clip attachment) slings
- Non-absorbent polyurethane padding




Specifications

Capacity	300kg
Graduations	100g
Dimensions	550mm (w) x variable height and length
Weight of scale	9.5kg

BODY COMPOSITION SCALES

Body Composition Scale 3552057




Features

- Measures Basal Metabolic rate, Total Body Water, Fat Free Mass, Fat Mass, fat %, weight and BMI
- Portable
- High Capacity
- Optional printer and carry case
- MDD approved

Specifications

Capacity	300kg
Graduations	100g
Dimensions	340mm x 450mm x 80mm
Weight of scale	9.5kg



Optional Carry Case



Body Composition Column Scale 3552058

Features

- Measures Basal Metabolic rate, Total Body Water, Fat Free Mass, Fat Mass, fat %, weight and BMI
- High capacity
- Large platform
- Optional digital or manual height measure
- Optional printer
- MDD approval

Specifications

Capacity	300kg
Graduations	100g
Dimensions	340mm x 450mm x 80mm
Weight of scale	12kg



Splashproof Swab Scale 3552077

Features

- Splashproof, IP65 rated
- Lightweight and compact design
- Over 1000 hours of continuous use
- Ideal for swab weighing
- Can be supplied with a flat plate - no scoop



Specifications

Capacity	6kg, 15kg, 30kg
Graduations	6kg: 1g / 15kg: 2g / 30kg: 5g
Dimensions	240mm (w) x 270mm (h) x 120mm (d)
Dimensions of plate	232mm x 179mm

CLASS
III
APPROVED

HOLD

TARE



Waterproof Stainless Steel Swab/ Organ Scale 3552080

Features

- Stainless steel scale
- High accuracy
- IP68 rated
- Red LED display
- Ideal for operating theatre (swab or organ weighing)
- Can be supplied with flat plate or a scoop fitted



Specifications

Capacity	6kg, 15kg, 30kg
Graduations	6kg: 1g / 15kg: 2g / 30kg: 5g
Dimensions	250mm (w) x 290mm (h) x 130mm (d)
Dimensions of plate	245mm (w) x 198.5mm (d)

CLASS
III
APPROVED

HOLD

TARE



Baby / Toddler Measuring Rod 3552067

Features

- Plastic mechanical measuring rod




Specifications

Increments	1mm
Dimensions	350mm x 800mm x 100mm
Measuring range	10-80cm
Weight	0.7kg

Digital Baby Height Measure 3552061

Features

- Fit 3552004 and 3552006 baby scales
- Aluminium digital height rod for measuring babies and toddlers
- Plastic baby tray
- Easy to assemble kit



Specifications

Increments	1mm
Dimensions	625mm x 355mm x 105mm
Measuring range	35-80cm
Weight	1.6kg

Portable Height Measure 3552070

Features

- Portable, light and easy to carry
- Moulded plastic design
- Six pieces, simple to assemble
- Complete with wall stabilisers
- Comes complete with robust corrugated plastic case



Specifications

Measuring range	0-2m x 1mm
Carry case dimensions	675mm x 365mm x 125mm

Telescopic Height Measure 3552064

Features

- Telescopic height measure made of anodised aluminium
- Can be fixed to the wall or mounted on column scales
- Incorporates fold down head piece for safety
- Patients height read at eye level



Specifications

Dimensions	35mm x 936mm x 26mm
Measuring range	60cm -2m
Increments	1mm

HEIGHT AND LENGTH MEASURES

Auto Digital Height Measure 3552071

Features

- Digital display
- Connects to Marsden column scales
- Marsden medical indicator (as on models 3552041 and 3552040)
- Offers automatic readout of weight, height and BMI on the indicator



Specifications

Dimensions	40mm x 30mm x 1780mm
Measuring range	1.2m - 2m
Increment	1mm
Weight	0.85kg
Power supply	1 x AAA Battery

Digital Telescopic Height Measure 3552065

Features

- Anodised aluminium digital height measure
- Attaches to Marsden column scales
- Data memory feature and data transfer capability



Specifications

Dimensions	40mm x 30mm x 1780mm
Measuring range	1.2m - 2m
Increment	1mm
Weight	0.85kg
Power Supply	1 x AAA Battery

Digital Wall-Mounted Height Measure 3552066

Features

- Anodised aluminium digital height measure
- Data memory feature and data transfer capability
- Suitable for wall mounting



Specifications

Dimensions	70mm x 80mm x 2100mm
Measuring range	1m - 2.1m
Increment	1mm
Weight	1.25kg
Power supply	1 x AAA Battery

Free Standing Height Measure 3552063

Features

- 2.54cm (1") square anodised aluminium height rod



Specifications

Base dimensions	350mm (w) x 50mm (h) x 380mm (d)
Minimum height	75cm
Maximum height	2.15m
Increment	1mm (1/8")

Body Measuring Rod 3552068

Features

- 2.54cm square anodised aluminium measuring rod
- Plastic slider fitted and foot support
- Available 1m or 2m



Specifications

3552068	1m
3552069	2m
Increment	1mm

Wall-Mounted Height Measure 3552067

Features

- 2.54cm square anodised aluminium height rod
- Wall fitting
- Attached to wall with two hinges (supplied)
- Plastic slider folds flat to the wall when not in use



Specifications

Minimum height	75cm
Maximum height	2.15m
Increment	1mm (1/8")

Baby Scale Rucksack 3552091

Features

- Suitable for 3552004, 3552006, 3552003 and 3552011
- Dimension 600mm (h) x 360mm (w) x 150mm (d)



Carry Cases

Baby Scale Carry Case 3552007

- Suitable for 3552004, 3552006, 3552011 and 3552003
- Dimensions (H x W x D): 590mm x 360mm x 140mm



Stand on Scales Carry Case 3552025

- Suitable for 3552016, 3552018, 3552020, 3552023 and 3552024
- Dimensions (H x W x D): 390mm x 350mm x 150mm



Carry Cases

Carry Case 3552029

- Suitable for 3552004, 3552027, 3552028
- Dimensions (H x L x W): 500mm x 650mm x 100mm



Hoist Weigher Carry Case 3552090 / 3552093



- Suitable for 3552056, 3552089
- 355209 Dimensions (H x L x W): 190mm x 610mm x 170mm
- 3552093 Dimensions (H x L x W): 190mm x 760mm x 170mm

Carry Case 3552095

Features

- Carry case suitable for hand-held baby scale
- Dimensions: 190mm (l) x 100mm diameter



Carry Case for Wheelchair Weigh Beams 3552051

Features


- Carry case suitable for wheelchair weigh beams
- Dimensions per beam (H x L x W):
190mm x 980mm x 120mm












Thermal Dot Line Printer 3552083

Features

- High resolution - 8 dots/mm
- Max. print width 48mm
- Print speed width 75mm/sec
- Print date, time, weight, height and BMI
- 3552103 has bluetooth connectivity



Product Specification Key

	Class Approval		Portable		Bluetooth
	Hold		Mains/ Non-Rechargeable		BMI Calculation
	Tare		Mains/ Rechargeable		Medical Devices Directive Approval



Active Mobility

Active Mobility Systems

Unit 8, 110-120 Silverwater Rd., SILVERWATER NSW 2128

P: 02 9649 2111 F: 02 9649 8506

W: www.activemobility.com.au E: sales@activemobility.com.au

