



Novacorr Hilo Health Bed



novacorr.com.au

Delivering the best quality and comfort in
health beds across Australia

Active Mobility | t. 02 9649 2111

novacorr

Hilo Health Bed



The Novacorr Hi-Lo Health Bed is designed to be raised or lowered at the touch of a button on the wireless hand control. Everyday tasks such as bathing, eating and taking medication also become much easier to manage for a carer or family member. When transferring from a wheel chair into/out of bed, the bed can be positioned at an ideal height to assist transferring. The Novacorr Hilo Health Bed has full hoist clearance.

Novacorr Healthcare can build a bed to suit the needs of everyone; weight, height and size is not an issue. We build our wireless remote control beds from small single to king and everything in between. The height of the bed can be made to an individual's height to assist with getting out of bed. The Adjustable Health bed can also be used purely as a luxury bed to make reading and watching TV more comfortable.

The back raise at the head of the bed can be raised from flat to ninety degrees (up right). This feature is useful if you suffer from the following conditions: back pain, hip/knee problems, breathing difficulties (including sleep apnoea, asthma, emphysema) GORD, arthritis. It also supports you into a sitting

position before getting out of bed. The leg raise at the foot of the bed can be raised to not only bend at the knee but keeps the whole leg elevated to support and improve circulation. This feature is useful if you suffer from swollen legs/ankles as it allows your legs to be higher than your heart. Using the leg raise feature can also help with hip/knee pain and will help alleviate back pain.

Novacorr Healthcare use state of the art Linak motors on our extensive range of beds. These motors are operated by using wireless remote controls. A simple touch of your control can deliver a perfect night's sleep and for your convenience no cumbersome wiring hanging from your bed. One of the greatest features on our beds is the ability to pre-set up to 3 different preferred positions. For example, by simply pressing the set button on your bed you can store your preferred sleeping position in the bed's computer memory. Likewise if you are advised by your Carer, Doctor or Occupational Therapist to keep your legs or back elevated to a predetermined height, that position can also be stored. Many of our clients also pre-set a position thus making it easier to exit your bed.

Mattresses (height of standard beds vary on mattress choice)

1. Inner Spring Pillow Top (10" or 254mm)



3. Memory Foam (8" or 203mm)



2. Contoured Latex (8" or 203mm)



4. Custom Made Mattress (prices may vary)

Features

- > Hilo Function to cover all Occupational Health and Safety needs
- > Head and leg elevation to allow the user optimal comfort
- > Steel frame design (custom made in our factory)
- > Weight tested up to 250kg (each side)
- > Heavy duty hospital grade motors
- > Multiple adjustable positions
- > Lockable hospital wheels
- > Foot and rail socks
- > Australian made and assembled
- > 10 year warranty

Sizes



Small Single
2'6" x 6'8" x 13"
762mm x 2032mm
x 330mm



Single
3'0" x 6'8" x 13"
914mm x 2032mm
x 330mm



King Single
3'6" x 6'8" x 13"
1066mm x 2032mm
x 330mm



Double
4'6" x 6'8" x 13"
1371mm x 2032mm
x 330mm



Queen
5'0" x 6'8" x 13"
1524mm x 2032mm
x 330mm



Dual Queen
2 x 2'6" x 6'8" x 13"
762mm x 2032mm
x 330mm



Dual King
2 x 3'0" x 6'8" x 13"
914mm x 2032mm
x 330mm



Dual Super King
2 x 3'6" x 6'8" x 13"
1066mm x 2032mm
x 330mm



Custom Made
Advise length, width &
height of custom bed
(prices may vary)

Optional Extras

1. Massage



2. Side Rails (may vary)



3. Overhead Poles



4. Side Support Rail



5. Valances



8. LCD Down Light



9. Hospital Wheels



6. Head and Foot Boards



7. Waterproof Antibacterial Breathable Mattress Protector



10. Straps



11. Bluetooth*



*To allow connectivity to hand control App for both Apple & Android

