



Novacorr Tilt and Lift Chair

novacorr.com.au

Delivering the best quality and comfort in
health chairs across Australia

Active Mobility | t. 02 9649 2111

Tilt and Lift Chair



Novacorr's range of custom made to measure recliner and sleeping chairs offer a level of comfort far beyond any chair purchased off the shelf at your local furniture store. Hand crafted by our highly skilled team in our own factory in Caloundra means we have the confidence in our chairs to offer a 10 year warranty on the frame, motors and mechanisms.

Novacorr's Tilt and Lift Chairs are custom built to the specifications of each individual client; we don't make chairs in standard small, medium or large sizes. This allows the client to have a chair that is not only comfortable, but suits their needs and requirements. Before any chair is made, a series of measurements are taken so your chair will fit you whilst ensuring there are no immediate pressure areas. Having your own custom made chair reduces the risk of back pain, shoulder strains and hip displacement whilst helping alleviate circulation problems and water retention plus many more benefits.

Once the frame on your chair has been made it is then covered in a selection of Dunlop memory foams. No two chairs are the same; the foam in the seating area will be made to your exact request. Your new chair is then covered in the fabric of your choice. Novacorr also supply a set of arm and neck covers and a full seat / leg overlay. These are designed to be removed and washed.

Novacorr's chairs are unlike any other chair, they are built to stay in a seated position as they recline. This allows for the client's feet to be elevated to a position where they are parallel with the heart. This will automatically help reduce swelling of the feet and legs whilst also boosting circulation. When the chairs are reclined they have a break in the knee, so when a client is lying down their knees are not locked at any point. A second motor can also be fitted which essentially turns the chair into a sleeping chair. The second motor enables the user to alter the back position without having to change the position of their legs.


Features

- > Marine ply chair frame
- > Linak and Dewert Motors being the leading electric actuator/motors manufacturers in the world
- > Dunlop commercial grade memory foams
- > Commercial grade webbing
- > Choice of Warwick, Windsor or Profile Fabrics
- > Emergency battery back-up
- > Commercial grade Dacron for extra comfort in the back rest
- > Removable arm covers, neck cushion and seat overlay helps extend the life of your chair
- > Side pockets for the convenience of placing items close
- > Reinforced steel mechanism

Standard choices

- > Fabric selection
- > Wing to be straight, half moon or elephant
- > Padding to be soft, medium or firm
- > Arm rest choices include skinny, standard or chubbo arms

We require certain measurements to begin a custom make of a chair. Please supply the following on purchase as per the diagram below:



- A - Floor to Seat
- B - Seat Width
- C - Seat Depth
- D - Back Height



Internal width
500mm



Seat depth
450 - 600mm



Lower leg
600mm



Vertical lift range
450mm - 650mm



Back recline
90 - 180 degrees



Foot elevation
to above heart



Weight capacity
130kg



Overall weight
up to 65kg




Warranty
10 years

Optional Extras

1. Massage



5. Extra Covers



2. Wheels



6. Waterproof Overlay



3. Hand Rail



7. Pressure Prevention Seating / Overlay



4. Latex Seating



8. Leather Upgrade



At Novacorr we consider ourselves to be problem solvers. If you have a request for anything specific to your chair let us know and we will do our best to make it work for you. Remember no one person is built the same and one size does not fit all!

