



Active Mobility

NDIS Registered Provider



## UNIVERSAL BOTTLE HOLDER

1669038

### FEATURES

- Adjusts to take most bottle sizes
- Easy to attach / adjust and re-position
- Swivel head for angle adjustment
- Adjustable mount can attach onto any tube / ledge up to 40mm thick
- Strong / secure grip
- Strong lightweight construction

### SPECIFICATIONS

#### Overall Dimensions

Width: 65mm  
Depth: 155mm  
Height: 90mm

#### Bottle Size

Ø55 - 70mm

#### Adjustable Grip

0 - 40mm (Will attach to round or flat surface)

#### Product Weight

160g

#### Maximum Load

1kg / 1 litre