


AERIS LITE

The AERIS Adjusta Lite manual wheelchair is a modern, adaptable wheelchair, suited to a large range of users.

Aeris is manufactured with aircraft grade Aluminium incorporating a host of features built into the design. AERIS is light weight, configurable & is fully crash tested to a user mass of 136kg Available in 16", 18", 20" seat widths



Colours



Calibre
Black

Sizes

16" Seat
18" Seat
20" Seat

Features



Aluminium frame

Aeris is manufactured with strong aircraft grade aluminium.



Weight Capacity

136 kg



Safety Tested

Crash tested and approved to 136kg



Configurable

Aeris incorporates a host of features built into the design.

Full Features

- Dual Brakes
- Solid 23" drive wheels. Pneumatic option is also available
- Solid 7.5" castors.
- Optional Accessories: Elevating legrest with calf support, Stump support and Headrest.
- Quick flip-up anti-tippers mechanism.
- Adjustable, Swing Away & Removable leg rests.
- Flip Back & Removable Armrests.
- Adjustable Centre of Gravity
- Height Adjustable Built-in push handles
- Footplates are angle adjustable.
- Removable or flip away, anti-tip wheels
- Swing away leg rests
- Adaptable backrest to suit other backrest options.
- Built-in backrest pocket

Specifications

| SPECIFICATIONS | AERIS LITE |
|-----------------------------|------------------------------|
| Weight Capacity | 136 kg |
| Seat Width | 40 cm, 46 cm, 49 cm |
| Armrest Height | 23 - 27 cm |
| Arm Pad Position | Fixed |
| Seat Height Adjustability | 44.5, 45, 47, 47.5, 50 cm |
| Seat Depth | 40, 42.5, 45, 47.5, 50 cm |
| Backrest Height | 40 - 46 cm (2cm increments) |
| Backrest Angle | 85° - 105° (5° increments) |
| Lower Leg Length | 40 - 48 cm (2 cm increments) |
| Footplate Length Adjustment | 19 -25 cm |
| Overall Height | 85 - 93 cm |
| Overall Length | 103 cm |
| Length excluding Footrest | 90 cm |
| Overall Width | is = to seat width +15 cm |
| Width Folded | 34 cm |
| Gross Weight | from 15 kg |
| Net Weight | from 9 kg |

Compliance

| COMPLIANCE | |
|-----------------------------|---|
| Australian Standards | Aeris complies with AS: 3695 / 3696 |
| ISO Crash Tested & Approved | Complies to ISO 7176-19:2008 as amended by EN121:2014 |

IMPORTANT : Always be sure to read your owner's manual throughly before operation.

Unit 8, 110-120 Silverwater Rd
Silverwater NSW 2128
Ph: 02 9649 2111

activemobility.com.au



AERIS
LITE