

# AERIS

## ADJUSTA




The AERIS Adjusta manual wheelchair series is a modern, adaptable wheelchair, suited to a large range of users.

The AERIS is manufactured with aircraft grade Aluminium incorporating a host of features built into the design. AERIS is light weight, highly configurable & is fully crash tested to a user mass of 136Kg. Available in 16", 18" 20" seat widths



### Colours



Calibre  
Black

### Sizes

16" Seat  
18" Seat  
20" Seat

### Features



**Aluminium frame**  
Aeris is manufactured with strong aircraft grade aluminium.



**Weight Capacity**  
136kg



**Safety Tested**  
Crash tested and approved to 136kg



**Highly Configurable**  
Aeris incorporates a host of features built into the design.

# Full Features

- Push to lock standard lock with optional drum brake kits
- Solid 23" drive wheels. Pneumatic option is also available
- Solid 7.5" castors.
- Accessories: elevating legrest with calf support, stump support and headrest.
- Removable and washable foam seating and backrest pad
- Quick flip-up anti-tippers mechanism.
- Adjustable Centre of Gravity
- Height Adjustable Built-in push handles
- Footplates are angle adjustable.
- Push to lock brakes
- Removable or flip away, anti-tip wheels
- Swing away leg rests
- Adaptable backrest to suit other backrest options.
- Built-in backrest pocket
- Adjustable, Swing Away & Removable leg rests.
- Flip Back & Removable Armrests.
- Tension Adjustable Backrest.

# Specifications

SPECIFICATIONS	AERIS ADJUSTA
Weight Capacity	136kg
Seat Width	40 cm, 46 cm, 49 cm
Footrest Extender	6 cm Footrest extender/joiner
Extendable Footrests	10.5 - 15.5 cm
Armrest Height	18 - 26 cm (2 cm increments)
Arm Pad Position	33 - 45 cm (4 cm increments)
Anterior Seat Height	42, 44.5, 45, 47, 47.5, 50 cm
Posterior Seat Height	41, 43, 45, 47, 47.5, 50 cm
Seat Depth	40, 42.5, 45, 47.5, 50 cm
Backrest Height	40 - 46 cm (2cm increments)
Backrest Angle	85° - 105° (5° increments)
Lower Leg Length	40 - 48 cm (2 cm increments)
Footplate Length Adjustment	10.5 -15.5 cm
Overall Height	85 - 93 cm
Overall Length	103 cm
Length excluding Footrest	90 cm
Overall Width	is = to seat width +15 cm
Width Folded	34 cm
Gross Weight	from 15.5 kg
Net Weight	from 9 kg

# Compliance

COMPLIANCE	
Australian Standards	Aeris complies with AS: 3695 / 3696
ISO Crash Tested & Approved	Complies to ISO 7176-19:2008 as amended by EN121:2014

**IMPORTANT** : Always be sure to read your owner's manual throughly before operation.

**Silverwater**  
 Unit 8, 110-120 Silverwater Rd  
 Silverwater NSW 2128  
 Ph: 9649 2111

**Newcastle**  
 7 Poynton Pl  
 Thornton, NSW 2322  
 Ph: 4047 3608