


AERIS ALIGNA

The **AERIS Aligna** manual wheelchair is a highly adjustable and ergonomic wheelchair that is suited to a range of users with its postural support and seated comfort.

The AERIS Aligna is crafted with a durable yet lightweight aluminum frame and waterproof upholstery, it offers reliability and ease of maintenance.



Colour Sizes Features



41cm
46cm
51cm



Aluminium frame
Manufactured with strong aircraft grade aluminium



Safety Tested
Crash tested and approved to 140kg



Weight Capacity
140kg



Configurable
Incorporates a host of features built into the design

Full Features

- **Padded Cushions:** optional pressure care upgrades available
- **Backrest:** gas- assisted tilt in space
- **Armrests:** adjustable
- **Leg rests:** elevating, swingaway and removable
- **Calf pads:** plastic flip up design
- **Footplates:** height and angle adjustable
- **Brakes:** push to lock wheel brakes and attendant operated drum brakes
- **Anti tip wheels**
- **Push handles**
- **Wheels:** quick release 61cm
- **Castors:** 200 x 25mm
- **Available in three widths:** 41cm, 46cm and 51cm

MOST FEATURES ARE ADJUSTABLE AND SOME REMOVEABLE

SPECIFICATIONS	
SEAT WIDTH	41/46/51cm
SEAT DEPTH	38-50cm
SEAT HEIGHT	41-52cm
BACK HEIGHT	50-65cm
ARM HEIGHT	18-29cm
PUSH HANDLE HEIGHT	98-116cm
OVERALL WIDTH	64/69/74cm
GROSS WEIGHT	38kg
NET WEIGHT	22kg
MAX USER WEIGHT	140kg

TGA Approved: 122644

Crash tested for safety

ISO 7176-19 :2008 / EN 121:2014

Unit 8

110-120 Silverwater Road

Silverwater NSW 2128

02 9649 2111

activemobility.com.au

sales@activemobility.com.au


Active Mobility

AERIS
ALIGNA