




FEATURES

Why Choose Eco Slings?

Eco Slings are a game-changing solution in patient handling—offering reliable support, minimised cross-contamination, and environmentally conscious disposal. Manufactured using 40% recycled polyester and engineered for strength, each sling supports a safe working load of up to 250kg. Designed with infection control in mind, Eco Slings combine safety with practicality.



Eco-Conscious Materials:

Made with 40% recycled polyester for reduced environmental impact



Robust & Reliable:

High safe working load of 250kg with durable stitching



Clear Patient Labelling:

Promotes efficient tracking and use and includes a wetness indicator



Universal Compatibility:

Fits most lifters and spreader bars



Size Options:

Available in Small, Medium, Large and Extra Large



Fixing Choices:

Offered in both Clip and Loop versions across all sizes



02 9649 2111

Unit 8

110-120 Silverwater Rd

Silverwater NSW 2128

activemobility.com.au