

# Hoist Attachments - General Purpose



## Standard Yoke

- 2 sling attachment points at extremity of curved yoke bar
- Use with general purpose (loop attachment) slings
- Non-absorbent polyurethane padding

**Available on:** Alto, Sonata & Sonata ELS

**Indications:** Suitable for the majority of patients

**Useful for:** All day-to-day patient lifting



Standard Yoke - SY02

## 4 Point Yoke

- 4 sling attachment points spread over a broad area
- Use with general purpose (loop attachment) slings
- Padded vinyl covers with velcro closure

**Available on:** Alto, Sonata & Sonata ELS

**Indications:** Bariatric clients, clients requiring a semi-supine lift to maintain upper body control

**Useful for:** Supine lifting devices eg. Jordan frame



Power Pivot Frame - PPF01

*Allegro*  
CONCEPTS PTY LTD

Innovation in Patient Handling



Active Mobility

## Hoist Attachments - General Purpose

The Allegro Concepts range of standard patient lifting hoists includes:

- Alto Lift (SWL 200kg)
- Sonata (SWL 150kg)
- Sonata ELS (SWL 150kg)

These three heavy duty hoists are available with a choice of 6 different sling attachment options.

Achieve quality patient transfers by selecting the fitting which best suits your client or healthcare facility.

Factors which may determine choice of sling attachment include:

- Physical attributes of patients eg. weight, disability
- Level of staff training eg. competence with complex lifting procedures
- Necessity to interface with existing equipment eg. beds, wheelchairs

### Tri-Hook Yoke

- 6 Sling attachment points at extremity of curved yoke bar
- Use with general purpose (loop attachment) slings
- Non-absorbent polyurethane padding

**Available on:** Alto, Sonata & Sonata ELS

**Indications:** Where the sling prescribed has more than the usual 4 attachment loops

**Useful for:** Facilities which have historically trained with this configuration



Tri Hook Yoke - SY03

### Hi-Lift Yoke


- 2 sling attachment points at extremity of slightly curved yoke bar
  - Use with general purpose (loop attachment) slings
  - Padded vinyl cover with velcro closure
- Available on:** Alto, Sonata & Sonata ELS
- Indications:** Where lifting height of specified hoist is insufficient (increases lifting height by 80mm)
- Useful for:** Reducing the loss of lifting height when installing an Allegro Universal Weighing Device



Hi-Lift Yoke - SY02HL

### Bariatric 4 Point Yoke - 3111040


- All day-to-day patient lifting
- Easily adaptable to any style of patient hoist
- Optimum patient space for larger patients
- 4 Sling attachment points
- Use with general purpose (loop attachment) slings
- Non-absorbent polyurethane padding
- 320kg SWL



Bariatric 4 Point Yoke - 3111040

### Bariatric Yoke Spreader Bar with Integrated Weigh Scale

- All day-to-day patient lifting
- Monitor patient wellbeing with no loss of hoist lifting height
- Easily adaptable to any style of patient hoist
- Optimum patient space for larger patients
- 4 Sling attachment points
- Use with general purpose (loop attachment) slings
- Non-absorbent polyurethane padding



Bariatric Yoke Spreader Bar with Integrated Weigh Scale - 3111039

## Active Mobility Systems Pty Ltd

Unit 8, 110- 120 Silverwater Road | Silverwater | NSW | 2128

☎ 02 9649 2111 | ✉ sales@activemobility.com.au

🌐 www.activemobility.com.au