

# SONATA 150

Compact Mobile Patient Lifting Hoist



*Allegro*  
CONCEPTS PTY LTD

Innovation in Patient Handling



Active Mobility

# SONATA 150

Sonata 150: A lightweight compact mobile patient lifting hoist with a safe working load of 150kg and an innovative and ergonomic design.

The Sonata 150 is an easy to use everyday lifter providing superior patient lift space through the use of SPACE LINK technology.



## SPACE LINK™

- breakthrough for optimized patient space
- enhances floor lift capacity

## Sonata 150 Specifications

Safe Working Load	150kg
Base Length	1005mm
External Base Width Closed	605mm
External Base Width Open	1070mm
Internal Base Width Closed	510mm
Internal Base Width Open	975mm
Base Width Adjustment	3 Position Manual
Under Bed Clearance	130mm
Minimum Height Yoke	420mm
Maximum Height Yoke	735mm
Average Lift Range*	870mm
Lifter Weight	41kg
Minimum Storage Height	1260mm
Under Base Height	55mm
Maximum Hoist Reach	875mm
Minimum Hoist Reach	555mm
Hoist Reach at Maximum Height	555mm
Mass of Heaviest Part	22kg (Mast assembly)

\* With general purpose large sling

The Allegro Concepts range of patient lifting equipment including lifters, slings and accessories has been designed and tested to comply with: AS ISO 10535-2011 Hoists for the transfer of disabled persons – Requirements and test methods.

## Active Mobility Systems Pty Ltd


Unit 8, 110- 120 Silverwater Road | Silverwater | NSW | 2128

☎ 02 9649 2111 | ✉ sales@activemobility.com.au

🌐 www.activemobility.com.au



Standard Yoke - 2211030



4 Point Yoke - 2211036



Standard Yoke with Integrated Weigh Scale - 2211037



Hi Lift Yoke - 2211032



Pivot Frame - 2211038



Quick Release Adapter - 2911039