



ATMOS

ULTRA LITE

As light as air, the Atmos ultra-lite, self-propel wheelchair weighs only 8kg and is one of the lightest wheelchairs in its class. Atmos combines the best in safety, portability and deluxe features in one ultra-lite wheelchair.



Colours



Sky Blue/
Black

Sizes

40cm Seat
46cm Seat
50cm Seat

Features



Ultra Light Weight
8kg¹



Weight Capacity
136kg



Safety Tested
Crash tested and
approved to 136kg



Ultra Portable
Folds in to a compact
unit suitable for most
car boots

Full Features

- Ultra-light weight frame
- Crash tested and approved to 136kg
- Reinforced cross brace frame design enhancing strength, durability, and performance.
- One-touch wheel release and collapsible frame folds down into one compact unit for travel.
- Removable/washable seat & backrest cushion
- Flip up footplates with integrated heel loop
- Cushioned seat distributes pressure, assists with posture and provides comfort
- Ergonomically designed comfort hand grips
- Large rear backrest pocket and armrest pocket for mobile phone
- Dual push to lock and user brakes
- Removable or flip away, anti-tip wheels
- Swing away leg rests

Specifications

SPECIFICATIONS	ATMOS SELF-PROPEL
Weight Capacity	136kg
Unit Weight ¹	8kg
Total Unit Width ²	10.5kg
Length	101cm
Height	92cm
Seat Width	40cm ,46cm 50cm
Seat Depth	40cm
Backrest Height	43cm
Folded Length	68cm
Folded Width	28cm
Folded Height	77cm
Height with Backrest Folded	73cm
Armrest Height to Seat	20cm
Height with Footplate Removed	76cm
Leg Rest Adjustment Range	10cm
Leg Rest Adjustment Range with Optional Extension	20cm
Self-Propel Wheel Size	22"
Front Castor Size	6" PU Filled
Rear Wheel Size	12" PU filled

Compliance

COMPLIANCE	
Australian Register of Therapeutic Goods (ARTG)	Atmos is entered in the ARTG 122644
ISO Crash Tested & Approved	Complies to ISO 7176-19:2008 as amended by EN121:2014

¹ Weight is calculated on a 46cm seat with leg rests, cushion and 22" wheel removed.

² Total unit weight is based on a 46cm seat width

^{*} Due to manufacturing tolerances and continual product improvement, all specifications can be subject to a variance of (+ or - 3%).

IMPORTANT : Always be sure to read your owner's manual thoroughly before operation.

Silverwater

Unit 8, 110-120 Silverwater Rd
Silverwater NSW 2128
Ph: 9649 2111

Newcastle

7 Poynton Pl
Thornton, NSW 2322
Ph: 4047 3608

activemobility.com.au



ATMOS