



ATMOS

ULTRA LITE

As light as air, the Atmos ultra-lite, self-propel wheelchair weighs only 8kg and is one of the lightest wheelchairs in its class. Atmos combines the best in safety, portability and deluxe features in one ultra-lite wheelchair.



Colours



Sky Blue/
Black

Sizes

40cm Seat
46cm Seat
50cm Seat

Features



Ultra Light Weight
8kg¹



Weight Capacity
136kg



Safety Tested
Crash tested and approved to 136kg



Ultra Portable
Folds in to a compact unit suitable for most car boots

Full Features

- Ultra-light weight frame
- Crash tested and approved to 136kg
- Reinforced cross brace frame design enhancing strength, durability, and performance.
- One-touch wheel release and collapsible frame folds down into one compact unit for travel.
- Removable/washable seat & backrest cushion
- Flip up footplates with integrated heel loop
- Cushioned seat distributes pressure, assists with posture and provides comfort
- Ergonomically designed comfort hand grips
- Large rear backrest pocket and armrest pocket for mobile phone
- Dual push to lock and user brakes
- Removable or flip away, anti-tip wheels
- Swing away leg rests

Specifications

| SPECIFICATIONS | ATMOS SELF-PROPEL |
|---|-------------------|
| Weight Capacity | 136kg |
| Unit Weight ¹ | 8kg |
| Total Unit Width ² | 10.5kg |
| Length | 101cm |
| Height | 92cm |
| Seat Width | 40cm, 46cm, 50cm |
| Seat Depth | 40cm |
| Backrest Height | 43cm |
| Folded Length | 68cm |
| Folded Width | 28cm |
| Folded Height | 77cm |
| Height with Backrest Folded | 73cm |
| Armrest Height to Seat | 20cm |
| Height with Footplate Removed | 76cm |
| Leg Rest Adjustment Range | 10cm |
| Leg Rest Adjustment Range with Optional Extension | 20cm |
| Self-Propel Wheel Size | 22" |
| Front Castor Size | 6" PU Filled |
| Rear Wheel Size | 12" PU filled |

Compliance

¹ Weight is calculated on a 46cm seat with leg rests, cushion and 22" wheel removed.

² Total unit weight is based on a 46cm seat width

^{*} Due to manufacturing tolerances and continual product improvement, all specifications can be subject to a variance of (+ or - 3%).

| COMPLIANCE | |
|---|---|
| Australian Register of Therapeutic Goods (ARTG) | Atmos is entered in the ARTG 122644 |
| ISO Crash Tested & Approved | Complies to ISO 7176-19:2008 as amended by EN121:2014 |

Warranty Information:

Frame: 2 Years

Upholstery: 12 Months

Other components (castors, brakes etc): 12 Months

IMPORTANT: Always be sure to read your owner's manual thoroughly before operation.

Unit 8, 110-120 Silverwater Rd
Silverwater NSW 2128
Ph: 02 9649 2111

activemobility.com.au



ATMOS