



# ATMOS

ULTRA LITE

As light as air, the Atmos ultra-lite transit wheelchair weighs only 7.5kg<sup>1</sup> and is one of the lightest wheelchairs in its class. Atmos combines the best in safety, portability and deluxe features in one ultra-lite wheelchair.



## Colours



Sky Blue/  
Black

## Sizes

40cm Seat  
46cm Seat  
50cm Seat

## Features



**Ultra Light Weight**  
7.5kg<sup>1</sup>



**Weight Capacity**  
136kg



**Safety Tested**  
Crash tested and approved to 136kg



**Ultra Portable**  
Folds in to a compact unit suitable for most car boots

# Full Features

- Ultra-light weight frame
- Crash tested and approved to 136kg
- Reinforced cross brace frame design enhancing strength, durability, and performance.
- One-touch collapsible frame folds down in to one compact unit for travel.
- Removable/washable seat & backrest cushion
- Flip Up Footplates with Integrated Heel Loop
- Cushioned seat distributes pressure, assists with posture and provides comfort
- Ergonomically designed comfort hand grips
- Large rear backrest pocket and armrest pocket for mobile phone
- Push to lock brakes
- Removable or flip away, anti-tip wheels
- Swing away leg rests

# Specifications

SPECIFICATIONS	ATMOS TRANSIT
Weight Capacity	136kg
Unit Weight <sup>1</sup>	7.5kg
Total Unit Weight <sup>2</sup>	9kg
Length	87cm
Height	92cm
Seat Width	40cm, 46cm & 50cm
Seat Depth	40cm
Backrest Height	43cm
Folded Length	65cm
Folded Width	28cm
Folded Height	73cm
Height with Backrest Folded	73cm
Armrest Height to Seat	20cm
Height with Footplate Removed	62cm
Leg Rest Adjustment Range	10cm
Leg Rest Adjustment Range with Optional Extension	20cm
Front Castor	6" PU Filled
Rear Wheel	12" PU filled

# Compliance

COMPLIANCE	
Australian Register of Therapeutic Goods (ARTG)	Atmos is entered in the ARTG 122644
ISO Crash Tested & Approved	Complies to ISO 7176-19:2008 as amended by EN121:2014

<sup>1</sup> Weight is calculated on a 46cm seat with leg rests & cushion removed.

<sup>2</sup> Total unit weight is based on a 46cm seat width

\* Due to manufacturing tolerances and continual product improvement, all specifications can be subject to a variance of (+ or - 3%).

**IMPORTANT :** Always be sure to read your owner's manual thoroughly before operation.

## Silverwater

Unit 8, 110-120 Silverwater Rd  
Silverwater NSW 2128  
Ph: 9649 2111

## Newcastle

7 Poynton Pl  
Thornton, NSW 2322  
Ph: 4047 3608