

CARBON

CRUISER



The Carbon Cruiser is a classic walker designed for those who value independence, style and smart functionality. Crafted from carbon fiber, it is incredibly lightweight yet durable, making everyday mobility smoother and more enjoyable.



Colours Available



Pheonix
Red



Obsidian
Black

Features



Portable
This walker is foldable and comes with a travel lock for easy transportation



Weight Capacity
150kg



Weight
5.5kg (without bag)

Full Features

- **Compact Side Fold Design:**
Easily folds sideways for a slim profile
- **Integrated Travel Lock:**
Keeps the walker securely folded during transport
- **Build-in Seat and Backrest:**
Comfortable seat and supportive back rest
- **Push to Lock Brakes:**
Simple, secure braking system for added safety and control
- **Soft Inverted Handles:**
Ergonomic grip to reduce strain and offer a natural hand position
- **Adjustable Hand Grips:**
Customisable height for a perfect fit and posture
- **Smooth Wheels:**
Glide effortlessly across indoor and outdoor surfaces
- **Included Bag:**
Carry your essentials with ease – stylish and practical

Main Technical Information	Carbon Cruiser
Material	Carbon Fiber
Max. Load	150kg
Max. Load of the Bag	6kg
Total Weight (without shopping bag)	5.5kg
Shopping Bag	0.6kg
Seat Width	460mm
Height of Hands Grips	780–980mm
Seat Height	580mm
Brakes	2 x maneuver brakes with parking function
Wheels (diameter x width)	Front: 200x36mm Rear: 200x36mm
Unfold Size	650x620x780mm
Fold Size	650x220x825mm
Package Size	660x230x880mm
Loading QTY	558PCS/HC
Total Weight with Carton Box	8.5kg
Standard Configuration	One shopping seat and one seat