

Azure  
DESIGNS

PORTABLE ROLLATOR WALKER

# ROVENA

**ROVENA** offers unmatched **comfort, stability, and convenience.**

Whether you're running daily errands or enjoying a leisurely stroll outdoors, its **durable aluminium frame** ensures **smooth** and **effortless** steering.



## Colours Available



Blue



Dark Red



Black

## Features



### Fold Design

The Rovena incorporates a compact front fold design that makes it easy to store and transport



**Weight Capacity**  
150kg



### Integrated Storage

Comes with a built-in storage bag conveniently tucked beneath the walker

# Full Features

- Sturdy steel frame, weighing under 9kg
- Durable, smooth-rolling 8-inch castors designed for both indoor and outdoor use
- Push-to-lock brakes for added safety — just push down to lock the wheels in place when stationary
- Thicker chair cushion ensures a more comfortable sitting experience, ideal for longer periods of use
- Thicker backrest for comfort & support
- Adjustable handle height from 840mm to 980mm, improving ergonomics and reducing strain
- Integrated baggage underneath the seat with a max load of 5kg
- Front fold-open mechanism allows for quick and compact storage
- Integrated fold-open chair provides access to storage and offers a reliable resting spot

| SPECIFICATIONS                      | ROVENA                                    |
|-------------------------------------|-------------------------------------------|
| Material                            | Aluminium                                 |
| Max. Load                           | 150kg                                     |
| Max. Load of the Bag                | 5kg                                       |
| Total Weight (without shopping bag) | 8.7kg                                     |
| Seat Width                          | 450mm                                     |
| Height of Hand Grips (adjustable)   | 840-980mm                                 |
| Seat Height                         | 530mm                                     |
| Brakes                              | 2 x maneuver brakes with parking function |
| Wheels (diamter x width)            | Front: 200x36mm<br>Rear: 200x36mm         |
| Unfold Size                         | 660x620x840mm                             |
| Fold Size                           | 910x620x220mm                             |

Unit 9  
110-120 Silverwater Rd  
Silverater NSW 2128  
02 9737 0144

*Azure*  
DESIGNS

azuredesigns.com.au  
sales@azuredesigns.com.au

**ROVENA**