



# CENTRO 2.0

Urban

By Active Mobility





# CENTRO 2.0

Terrain

# Introducing the NEW CENTRO 2.0



## Better than ever...

Proudly designed and manufactured in Australia, the Centro 2.0 is here.

The new Centro 2.0 has many modern updates without losing any of its robust and reliable design features.

## FEATURES

- Available in Adult, Paediatric, Bariatric Version and XT
- 50° seat tilt degree sensor
- New shear reduction back recline
- Memory Positioning
- Multiple wheel rim and tyre options
- Tuned APT suspension
- 80mm ground clearance
- Electronic Brakes with manual park
- Vehicle tie down points
- Crash Tested
- Extensive colour range to choose from





Low-Shear Recline



## Low Shear Recline

Optional Low Shear recline available on standard and vertical lift models



Tilt

## Posterior Tilt

For comfort and convenience the standard seat has 50° posterior tilt.

The tilt angle sensor indicator will display on the control screen in degrees as angle is adjusted.




## Vertical Lift

Users of the Centro 2.0 vertical lift can raise their seat by 30cm, perfect for social settings



Raise

# FEATURES



### APT Suspension

6-wheel tuned suspension allows for the smooth transverse of kerbs.



### Customised Seating

The Centro 2.0 can adapt to any seating system.




### Wheel Choices

A wide range of wheel options such as wide terrain tyres. (see next page)



### Colour Choices

Choose from an extensive colour palette or a design of your choice.



### Memory Positioning

The Centro 2.0 can be programmed to hold up to 15 memory positions.



### Elevating Legrests

Power elevating, centre mount and swing away legrests





**Tackle Outdoor Obstacles with ease.**


# CENTRO 2.0

## Terrain Wheels

### Terrain

Best suited for the more adventurous indoor & outdoor user.

Wider footprint, Higher Tyre profile and more aggressive tread pattern provide superior Traction & outdoor performance, with minimal impact on indoor use.



### Extreme Terrain (XT)

Full Outdoor Knobby Low-pressure Tyres fitted to alloy Wide mag style wheels on the Centro XT model.

The ultimate outdoor tyre & wheel combo that gives Maximum Traction and shock absorbing characteristics.





**Your wheels, your choice**

# CENTRO 2.0

## Urban Wheels

Urban and Terrain wheels can be easily interchanged giving the best of both worlds



Standard Grey



Patterned Cover



Polished Alloy



The urban wheels can be used as general purpose for both indoor & outdoor activities.

These wheels provide the smallest footprint, maximum manoeuvrability with balanced performance and traction for everyday use.

Available as Pneumatic or Puncture proof PU filled.

### Castors

7" 8" & 9" diameter depending on Drive wheel configuration available in Polished mag style or Black





**Explore Colour & Design like never before..**

Leading technology allows for a an unlimited range of colour choices and design patterns, allowing you to keep your personal sense of style!



## Custom Designs



# CENTRO 2.0 Specifications

SPECIFICATIONS	Centro 2.0 Urban & Terrain
Recommended User Weight	Up to 175 kg (as per AS 3935)- full compliance
Weight	85kg (base)
Overall Width	620mm (wheel base)
Overall Length	900mm (with footplates removed)
Seat Widths	31 - 58 cm (custom sizes available)
Seat Depths	31 - 58 cm (custom sizes available)
Seat Height	46 cm (top of rail using 0-46 set up)
Low Seat Height Option	42 cm
Push Handle Height Range	100 - 105 cm top of push handle (custom heights available)
Back Height	1030 mm top of push handle (custom height available)
Seat Tilt	4 - 50°
Arm Rests	Flip up (adjustable)
Leg Rests	Power elevating and swing away detachable.
Frame Colour	Wide choice
Brakes	Electronic and manual park
Rear Castors	7" pneumatic (options: 6" or 8")
Front Castors	8" pneumatic (option: 7")/ or 8" solid fill
Drive Wheels	Urban: 3.00-8 (14")- options: 14" or 12" wide Terrain: 4.00-6 Pneumatic
Controls	LinX (option: R-net)
Battery	70-90 amp/hr, Sealed
Range	40 km (max- with 70 amp/hr)
Crash Tested	YES ( AS 3695, ISO 7176-19:2008, AS/NZS 3696.19:2009)

**Book a Trial online :** [activemobility.com.au](http://activemobility.com.au)

Or contact [sales@activemobility.com.au](mailto:sales@activemobility.com.au)



02 9649 2111



Unit 8, 110-120 Silverwater Road,  
Silverwater, NSW, 2128