



Active Mobility

Configura[®] Bariatric[™]





Configura® Bariatric brings Configura® technology to bariatric seating. This offers unprecedented adjustability, flexibility and peace of mind for the prescriber.

It is also available with a full width alternating air cushion to meet the needs of users at high risk of pressure damage.

EASY ADJUSTMENT before and after delivery to suit varied and changing needs

QUICK and easy to recycle via removable, washable covers

AFFORDABLE, providing more consistent budget management and significant savings on special chairs

MODEL	Heights	Widths	Depths
250kg	16", 18"	24", 26", 28"	18", 20", 22"
300kg	18"	30", 32", 34"	18", 20", 22"



Dual Motor

Configura® Bariatric's Dual Motor operation allows independent adjustment of the legrest elevation and backrest recline.



Tilt-in-space (with independent legrest)

Tilt-in-space improves pelvic stability and posture, pressure distribution, circulation, makes hoisting easier and reduces the tendency to slip forward in the chair.



3

Width adjustment

The Configura® Bariatric range makes selecting the chair width easy and trouble-free via removable armrests. This allows three different widths. to be achieved within the same chair giving you peace of mind and allows adjustment for changing needs after delivery.



4

Supportive armrests for easy transfers.



5

100kg weight capacity on legrest.



6

Breathable, wipeable fabric on seat base and backrest.

8



Pressure cushion options

A wide range of pressure cushions to suit different user needs:



Castellated



Visco



Cushionair cushion used with battery powered pump

7

Posture options



Multi-adjustable Pillow Backrest

Each pillow is height adjustable and filling can be customised to the user. Offering mild to moderate postural support, the Multiadjustable Pillow Backrest is highly effective for users with a kyphosis or a scoliosis. Excellent for postural accommodation. Supplied as standard with VP fabric for comfort and pressure reduction.



2 x Lateral Supports

Lateral Supports to fit underneath the backrest pillows to provide lateral support combined with the comfort and adjustability of the Multiadjustable Pillow Backrest.



Profiled Headrest

Suitable for users with reduced head control. Fits over the shoulders.



The Configura Concept

Assess

Using the Configura® Assessment Form, select the correct chair including options for pressure relief and postural support.

Deliver

Chair is delivered within 2 working days, configured according to specification chosen on Configura® Assessment Form. Adjustments can be made on-site to ensure maximum comfort and effectiveness.

Adjust

If the user's needs change over time, adjustments can be made such as increasing pressure care or providing more lateral support. In some cases, this will remove the need to issue a new chair.

Recycle

Configura® can be easily and cost effectively recycled as all contact areas can be removed for cleaning or replacement. The chair can then be re-configured for the next user.....

Training

We are committed to providing training and support to make you competent and confident in issuing Configura chairs.

There are two ways we can help prescribers:

Product awareness days - we visit your organisation to carry out detailed product training sessions covering the functions and issuing of the chair.

Professional seminars - day long seminars delivered by professional Occupational Therapists to raise awareness of seating issues such as pressure and posture management.



At Active Mobility we source quality, affordable and innovative products globally and have them available for fast delivery. But we do more than just provide products.

We give the knowledge and support to help our customers grow their business and deliver better life care.



Active Mobility

ACTIVE MOBILITY SYSTEMS PTY LTD

Unit 8/110-120 Silverwater Road

Silverwater NSW 2128

Ph: 02 9649 2111 | @: sales@activemobility.com.au

www.activemobility.com.au