




Colibri

Stylish, simple and micro lightweight. The ultra portable Colibri can be disassembled in under a minute for transportation and is designed for those who enjoy a totally independent lifestyle.



Colours Available



Features



Seat
Mid Back Soft Padded Seat



Wheels
4 Wheel



Weight Capacity
136kg



Lighting
Front and Side Reflectors



Portable
Dismantle/ Reassemble

Full Features

- Effortless Dismantling in Less Than One Minute Without the Need for Tools
- Easy Fit Into Car Boots
- Compact Design to Maneuver in Tight Spaces
- Charge Battery Directly on Scooter or Detached From It
- Optional 'Easy Fit Seat' Suspension Kit
- Unlocking Lever for Swiveling and Removing Seat
- Soft Padded Seat for Comfort
- Ergonomic Handle Bar
- Fold Arm Rests
- Light Weight and Portable
- Front and Side Reflectors
- Easy Operating Console
- Front Basket
- Horn

Specifications

SPECIFICATIONS	INVACARE COLIBRI
Backrest Width	46.5cm
Seat Depth	40cm
Floor to Seat Height	Minimum: 37cm, Maximum: 42cm
Backrest Height	38.5cm
Driving Unit Width	50.5cm
Total Length	101cm
Empty Weight	49kg
Max User Weight	136kg
Battery Capacity	2 x 18 Ah
Turning Radius	220cm
Max Climbable Obstacle Height	4.5cm
Ground Clearance	5cm
Maximum Safe Slope	6° Degrees / 10° Degrees
Maximum Driving Range ¹	16km
Maximum Driving Speed	8km/h
Motor	200 W

¹Driving range varies with the factors shown below:

- (a) The weight of the occupant
- (b) Level ground & ambient temperature
- (c) Battery condition / Type of charger
- (d) Driving Style

Warranty Information :

Scooters: 5 year warranty
 Batteries and Chargers: 12 Months
 Scooter Accessories: 12 months on all accessories

IMPORTANT : Always be sure to read your owner's manual thoroughly before operating your machine.

WARNING : Do not exceed 10km/ph on public roads or paths.