



Comet

Alpine +

Robust and powerful, the Comet Alpine+ gives you the freedom to discover the great outdoors with minimal effort and in continuous comfort and safety.



Colours Available



Jade Green



Allum Blue



Electric Blue



Garnet Red



Carmine Red



Crystal Copper



Luxor Beige



Ionic Chrome



Diamond Silver



Marble White

Features



Seat

High Back
Contoured Captain
Seat with Electronic
Seat Riser



Wheels

4 Wheel



Weight Capacity

160kg
Bariatric Suitable



Lighting

LED Headlights,
Indicators, Hazard Lights

Full Features

- Ability to Drive Across Various Terrains
- Powerful and Robust
- Superior Suspension System
- Electronic Seat Riser to Achieve Preferred Driving Height
- In-Built Splash Guards to Protect From Water and Dirt
- Directional Stability
- Automatic Speed Reduction When Driving Around Corners
- Clear, Easy to Read LCD Dashboard
- High Back Contoured Seat with Headrest
- Front Basket
- LED Headlights
- Hazard Lights / Indicators
- Excellent Outdoor Performance

Specifications

SPECIFICATIONS	INVACARE COMET ALPINE+
Backrest Width	51cm
Seat Depth	46cm
Floor to Seat Height	Minimum: 44cm, Maximum: 51cm
Backrest Height	63cm
Driving Unit Width	66.5cm - 68.5cm
Total Length	145cm
Empty Weight	143kg
Max User Weight	160kg
Battery Capacity	2 x 12-volt / 75Ah
Turning Radius	275cm
Max Climbable Obstacle Height	10cm
Ground Clearance	10cm
Maximum Safe Slope	10° Degrees
Maximum Driving Range ¹	58km
Maximum Driving Speed	10km/h
Motor	650 W / 1800 W

¹Driving range varies with the factors shown below:

- (a) The weight of the occupant
- (b) Level ground & ambient temperature
- (c) Battery condition / Type of charger
- (d) Driving Style

Warranty Information :

Scooters: 5 year warranty

Batteries and Chargers: 12 Months

Scooter Accessories: 12 months on all accessories

IMPORTANT : Always be sure to read your owner's manual thoroughly before operating your machine.

WARNING : Do not exceed 10km/ph on public roads or paths.