



Comet

Ultra

Ideal for those with larger body shapes, the Comet Ultra offers a secure, supportive and powerful driving experience without compromising on style or functionality.



Colours Available



Jade Green



Allum Blue



Electric Blue



Garnet Red



Carmine Red



Crystal Copper



Luxor Beige



Ionic Chrome



Diamond Silver



Marble White

Features



Seat

High back contoured ultra seat with lumbar support



Wheels

4 Wheel



Weight Capacity

220kg
Bariatric suitable



Lighting

LED headlights, indicators, hazard lights

Full Features

- Ability to drive across various terrains
- Powerful and robust for optimum safety and stability
- Special contoured ultra seat developed by therapists to ensure optimum seating position
- Reinforced height & width adjustable armrests
- Superior suspension to offer added comfort for larger individuals
- Reinforced seat base
- Clear, easy to read LCD dashboard
- Front basket
- LED headlights
- Hazard lights / indicators
- Designed for those with larger body shapes

Specifications

SPECIFICATIONS	INVACARE COMET ULTRA
Maximum Weight Capacity	220kg
Maximum Speed	10km/h
Maximum Range ¹	49km
Length	145cm
Width	66.5cm - 68.5cm
Scooter Weight	148kg
Ground Clearance	10cm
Seat to Ground	Minimum: 44cm, Maximum: 51cm
Turning Radius	275cm
Battery Requirements	2 x 12-volt / 75Ah
Seating	Contoured Ultra Seat, Depth 53.5cm
Motor	650 W / 1800 W

¹Driving range varies with the factors shown below:

- (a) The weight of the occupant
- (b) Level ground & ambient temperature
- (c) Battery condition / Type of charger
- (d) Driving Style

Warranty Information :

Scooters: 5 year warranty

Batteries and Chargers: 12 Months

Scooter Accessories: 12 months on all accessories

IMPORTANT : Always be sure to read your owner's manual thoroughly before operating your machine.

WARNING : Do not exceed 10km/ph on public roads or paths.

Silverwater

Unit 8, 110-120 Silverwater Rd
Silverwater NSW 2128
Ph: 9649 2111

Newcastle

7 Poynton Pl
Thornton, NSW 2322
Ph: 4047 3608

activemobility.com.au

