




Comet Ultra

Ideal for those with larger body shapes, the Comet Ultra offers a secure, supportive and powerful driving experience without compromising on style or functionality.



Colours Available



Features



Seat
High back contoured ultra seat with lumbar support



Weight Capacity
220kg
Bariatric suitable



Wheels
4 Wheel



Lighting
LED headlights, indicators, hazard lights

Full Features

- Ability to drive across various terrains
- Powerful and robust for optimum safety and stability
- Special contoured ultra seat developed by therapists to ensure optimum seating position
- Reinforced height & width adjustable armrests
- Superior suspension to offer added comfort for larger individuals
- Reinforced seat base
- Clear, easy to read LCD dashboard
- Front basket
- LED headlights
- Hazard lights / indicators
- Designed for those with larger body shapes

Specifications

SPECIFICATIONS	INVACARE COMET ULTRA
Maximum Weight Capacity	220kg
Maximum Speed	10km/h
Maximum Range ¹	49km
Length	145cm
Width	66.5cm - 68.5cm
Scooter Weight	148kg
Ground Clearance	10cm
Seat to Ground	Minimum: 44cm, Maximum: 51cm
Turning Radius	275cm
Battery Requirements	2 x 12-volt / 75Ah
Seating	Contoured Ultra Seat, Depth 53.5cm
Motor	650 W / 1800 W

¹Driving range varies with the factors shown below:

- (a) The weight of the occupant
- (b) Level ground & ambient temperature
- (c) Battery condition / Type of charger
- (d) Driving Style

Warranty Information :

Scooters: 5 year warranty
 Batteries and Chargers: 12 Months
 Scooter Accessories: 12 months on all accessories

IMPORTANT : Always be sure to read your owner's manual thoroughly before operating your machine.

WARNING : Do not exceed 10km/ph on public roads or paths.

Silverwater
 Unit 8, 110-120 Silverwater Rd
 Silverwater NSW 2128
 Ph: 9649 2111

Newcastle
 7 Poynton Pl
 Thornton, NSW 2322
 Ph: 4047 3608

activemobility.com.au

