



# Pegasus

## Pro

Comfort, Safety and Reliability with added performance, the Pegasus Pro is the ideal choice for users who demand more.



## Colours Available



Jade Green



Allum Blue



Electric Blue



Garnet Red



Carmine Red



Crystal Copper



Luxor Beige



Ionic Chrome



Diamond Silver



Marble White

## Features



**Seat**  
High Back  
Contoured Captain  
Seat with Head Rest



**Wheels**  
4 Wheel



**Weight Capacity**  
160kg  
Bariatric Suitable



**Lighting**  
LED Headlights,  
Indicators, Hazard Lights,  
Brake Lights

# Full Features

- Advanced Suspension System
- Smooth Ride Over All Terrains
- Powerful Motor
- Pneumatic Tyres
- Additional Electronics and Motor Protection to Guard Against Water and Corrosion
- Easy Adjustable Tiller
- High Back Contoured Seat with Headrest
- On-Board Electronic Speed Reduction System Reduces Speed When Driving Around Corners
- Hand Brakes for Immediate Braking
- Clear, Bright Brake Lighting
- Easy to Read LCD Dashboard
- Front Basket
- LED Headlights
- Hazard Lights / Indicators

# Specifications

SPECIFICATIONS	INVACARE PEGASUS PRO
Backrest Width	51cm
Seat Depth	47cm
Floor to Seat Height	Minimum: 44cm, Maximum: 51cm
Backrest Height	49cm
Driving Unit Width	65cm
Total Length	132cm
Empty Weight	136kg
Max User Weight	160kg
Battery Capacity	2 x 12-volt / 75Ah
Turning Radius	275cm
Max Climbable Obstacle Height	10cm
Ground Clearance	10cm
Maximum Safe Slope	10° Degrees
Maximum Driving Range <sup>1</sup>	52km
Maximum Driving Speed	10km/h
Motor	550 W / 1300 W

<sup>1</sup>Driving range varies with the factors shown below:

- (a) The weight of the occupant
- (b) Level ground & ambient temperature
- (c) Battery condition / Type of charger
- (d) Driving Style

## Warranty Information:

Scooters: 5 year warranty

Batteries and Chargers: 12 Months

Scooter Accessories: 12 months on all accessories

**IMPORTANT :** Always be sure to read your owner's manual thoroughly before operating your machine.

**WARNING :** Do not exceed 10km/ph on public roads or paths.