




Pegasus Pro

Comfort, Safety and Reliability with added performance, the Pegasus Pro is the ideal choice for users who demand more.



Colours Available



Features



Seat
High Back
Contoured Captain
Seat with Head Rest



Weight Capacity
160kg
Bariatric Suitable



Wheels
4 Wheel



Lighting
LED Headlights,
Indicators, Hazard Lights,
Brake Lights

Full Features

- Advanced Suspension System
- Smooth Ride Over All Terrains
- Powerful Motor
- Pneumatic Tyres
- Additional Electronics and Motor Protection to Guard Against Water and Corrosion
- Easy Adjustable Tiller
- High Back Contoured Seat with Headrest
- On-Board Electronic Speed Reduction System Reduces Speed When Driving Around Corners
- Hand Brakes for Immediate Braking
- Clear, Bright Brake Lighting
- Easy to Read LCD Dashboard
- Front Basket
- LED Headlights
- Hazard Lights / Indicators

Specifications

SPECIFICATIONS	INVACARE PEGASUS PRO
Backrest Width	51cm
Seat Depth	47cm
Floor to Seat Height	Minimum: 44cm, Maximum: 51cm
Backrest Height	49cm
Driving Unit Width	65cm
Total Length	132cm
Empty Weight	136kg
Max User Weight	160kg
Battery Capacity	2 x 12-volt / 75Ah
Turning Radius	275cm
Max Climbable Obstacle Height	10cm
Ground Clearance	10cm
Maximum Safe Slope	10° Degrees
Maximum Driving Range ¹	52km
Maximum Driving Speed	10km/h
Motor	550 W / 1300 W

¹Driving range varies with the factors shown below:

- (a) The weight of the occupant
- (b) Level ground & ambient temperature
- (c) Battery condition / Type of charger
- (d) Driving Style

Warranty Information:

Scooters: 5 year warranty

Batteries and Chargers: 12 Months

Scooter Accessories: 12 months on all accessories

IMPORTANT : Always be sure to read your owner's manual thoroughly before operating your machine.

WARNING : Do not exceed 10km/ph on public roads or paths.