

**NEW**

  
**Active Mobility**



# TDX Combo

Exclusive to Active Mobility, the Invacare TDX Combo sets a new standard in mid-wheel drive performance, combining advanced stability, intelligent control and fully integrated seating to support complex clinical needs.



## Colour



Wet Black

## Features



**Top Speed**  
10km/h



**SWL**  
136kg



**Legrest Options**  
Centre-Mount or  
Swing-Away



**Weight Capacity**  
136kg



**Tilt-in-Space**  
50°



**Seating Included**  
Headrest, seat belt,  
cushion and backrest

Exclusive to Active Mobility, the Invacare TDX Combo is designed to move with you. This high-performance mid-wheel drive power wheelchair offers smooth, confident driving, easy manoeuvrability and all-day comfort. With advanced Ultra Low Matrx seating, including a seat belt, back support and headrest, and intuitive LiNX controls, it adapts to your body and lifestyle so you can feel supported, comfortable and in control wherever your day takes you.

## Key Features:

- Tilt-in-space - 50°
- LiNX Control System
- Retractable Joystick Control Mount
- Swing Away Legrests - 70°
- Centre-Mount Legrest
- LCD LiNX Touch Screen Joystick Module
- 50Ah Gel Batteries
- Battery Charger Included
- 14" x 3" Black Tyres with Gel Foam Inserts
- LED Lights and Indicators
- Wheelchair Transport Brackets
- Seat Belt
- Matrx E2 Backrest
- Matrx Libra Cushion
- Invacare Matrx Elan Headrest
- Crash Tested for Safety

Specifications	
Seat Width	40.5cm - 61cm
Seat Depth	40.5cm - 58cm
Seat Height	43.82cm
Base Length	90.17cm
Base Width	60.96cm, 64.77cm
Turning Radius	20°
Max User Weight	136kg
Total Weight	From 136kg
Max Safe Slope	9°
Motor Type	4 Pole
Top Speed	10km/h
Compliance	AS/NZS 3695.2   ISO 7176-19
Shroud Colour	Wet Black
Warranty	Frame 5 years Components 1 year
TGA Approved	122644