



Jovy II

Ultra-Light Carbon Fibre Wheelchair

The ultimate in portable lightweight folding Power Wheelchairs.

This revolutionary design removes the requirement for a separate battery unit, further streamlining the chair design and reducing its total weight.



Colours



Size

49cm Seat

Features



Ultra Light Weight
21kg¹



Weight Capacity
120kg



Driving Range
high performance motors that can take you further for longer



Ultra Portable
One touch folding capability

Full Features

- Ultra-light weight frame
- Ultimate outdoor performance and a greater driving range as a result of high performance motors that can take you further for longer
- One touch folding capability and a battery unit that fits directly on the chair frame
- Generously wide and cushioned seat distributes pressure, assists with posture and provides comfort
- Ultra-light aircraft grade carbon fibre construction
- Removes the requirement for a separate battery unit, further streamlining the chair design and reducing its total weight.

Specifications

SPECIFICATIONS	JOVY II CarbonFibre
Weight Capacity	120kg
Unit Weight ¹	21kg
Dimensions Unfolded (L x W x H)	103 x 62 x 94cm
Dimensions Folded (L x W x H)	62 x 40 x 87cm
Seat Width	49cm
Seat Depth	40cm
Backrest Height	50cm
Turning Radius	85cm
Drive Model	Rear wheel drive
Motor	2x 180W brushless motor
Battery	2x 24V 6Ah lithium battery
Max Speed	6km/h
Drive Range	24km
Front Wheel Size	8" PU solid wheels
Rear Wheel Size	12.5 x 2.25" PU solid wheels
Climbing Slope	6°
Charging Time	6-8 hours



¹ Weight is calculated on a 49cm seat with battery units removed.
² Total unit weight is based on a 49cm seat width
^{*} Due to manufacturing tolerances and continual product improvement, all specifications can be subject to a variance of (+ or - 3%).

IMPORTANT : Always be sure to read your owner's manual thoroughly before operation.

Silverwater
 Unit 8, 110-120 Silverwater Rd
 Silverwater NSW 2128
 Ph: 9649 2111

Newcastle
 7 Poynton Pl
 Thornton, NSW 2322
 Ph: 4047 3608