

Kingsdene Change Table - **Adult**



The Kingsdene is designed to address a wide range of manual handling issues to benefit both patient and carers.

SKU: 3202005
ARTG: 312129



Colours




Blue

Size

Length
1970mm
Width
710mm

Features



Cot Sides
Fitted upholstered cot sides



Adjustable
Utilises a high grade electric actuator



Weight Capacity
200kg



Twin Lift
Strong twin-lift mechanism ensures rigidity and stability during treatment

Full Features

- Designed to address a wide range of manual handling issues to benefit both patient and carers.
- The height adjustment allows for comfortable handling of the patient to the correct positional height during procedures with enhanced safety and care.
- Fitted upholstered cot sides, where both sides are collapsible.
- Utilises a high-grade electric actuator that's driven to enable full adjusted product hi-lo.
- Strong twin-lift mechanism ensures rigidity and stability during treatment.
- Four 5" castors, all with brakes.
- **Complies to ISO 9001:2015, ISO13485:2016, ISO14001:2015, ISO 45001:2018 CE**

Specifications

Specifications	
Construction	High quality steel sections, and sheet construction utilizing the latest manufacturing techniques throughout.
Upholstery	Manufactured in line with using high quality foam, lining with anti-bacterial and fire retardant vinyl
Overall Width	710 mm
Overall Length	1970 mm
Cot Sides	Upholstered and Collapsible
Height	525- 975 mm
Coating	Expoxy powder coating
Colour	Blue
Item Weight	Approx. 63 kg

Examples



Silverwater

Unit 8, 110-120 Silverwater Rd
Silverwater NSW 2128
Ph: 9649 2111

IMPORTANT : Always be sure to read your owner's manual and installation guide thoroughly before operation.

activemobility.com.au


Active Mobility


ARCHABILITY