

Kingsdene Change Table - Junior



The Kingsdene is designed to address a wide range of manual handling issues to benefit both patient and carers.

SKU: 3202006
ARTG: 312129



Colours




Blue

Size

Length
1460mm
Width
710mm

Features



Cot Sides
Fitted upholstered cot sides



Adjustable
Utilises a high grade electric actuator



Weight Capacity
200kg



Twin Lift
Strong twin-lift mechanism ensures rigidity and stability during treatment

Full Features

- Designed to address a wide range of manual handling issues to benefit both patient and carers.
- The height adjustment allows for comfortable handling of the patient to the correct positional height during procedures with enhanced safety and care.
- Fitted upholstered cot sides, where both sides are collapsible.
- Utilises a high-grade electric actuator that's driven to enable full adjusted product hi-lo.
- Strong twin-lift mechanism ensures rigidity and stability during treatment.
- Four 5" castors, all with brakes.
- **Complies to ISO 9001:2015, ISO13485:2016, ISO14001:2015, ISO 45001:2018 CE**

Specifications

| Specifications | |
|-----------------------|---|
| Construction | High quality steel sections, and sheet construction utilizing the latest manufacturing techniques throughout. |
| Upholstery | Manufactured in line with using high quality foam, lining with anti-bacterial and fire retardant vinyl |
| Overall Width | 710 mm |
| Overall Length | 1460 mm |
| Cot Sides | Upholstered and Collapsible |
| Height | 525- 975 mm |
| Coating | Expoxy powder coating |
| Colour | Blue |
| Item Weight | Approx. 63 kg |

Examples



Silverwater

Unit 8, 110-120 Silverwater Rd
Silverwater NSW 2128
Ph: 9649 2111

IMPORTANT : Always be sure to read your owner's manual and installation guide thoroughly before operation.

activemobility.com.au



Active Mobility



ARCHABILITY