

SKU:  
2358000

# Monarch Portable Hoist

The Monarch Portable Ceiling Hoist can be used as a flexible alternative to fixed overhead hoists.



## Accessories\*

Floor Trolley: 2322512

Flexi-Link : 2322503

\*purchase separately

## Features



**Secure**  
Spreader Bar  
with secure sling  
mechanism



**Long Lasting  
Battery**  
Increased battery  
capacity



**Weight Capacity**  
200kg




**Portable**  
Ideal for  
transporting with  
the user

## Full Features

- + Ultra-Lightweight, weighing only 4.8kg
- + Comes standard with a unique spreader bar which includes a secure sling mechanism
- + Perfect for rental properties, multi-room, temporary or short-term use
- + Portable, which makes it ideal for transporting with the user
- + Easy Transport Handle
- + Available with telescopic reacher
- + Is able to lift the industry standard of 200kg.

## Specifications

- + Max Capacity : 200kg
- + Function : Patient Lifting, Room to Room Transfer, Temporary Use
- + Weight : 4.8kg
- + Battery Operation : 36 V lead acid batteries
- + Lifts per Charge : 30 to 40 cycles
- + Lifting Speed : 4.5cm/second
- + **Soft Start/Stop : YES**
- + **Emergency Stop/Lower : YES**
- + **Track Options** : H System, Turntable, Gate, Curves, Room-to-Room, Freestanding Gantry
- + **Track Compatibility** : Compatible with all Monarch tracking and other leading track manufacturers



## Applications

The Monarch Portable Hoist is Compatible across a range of applications, these include:



Nursing Homes



Hospitals



In-Home Care



Hotels/Rental Properties



Schools




Leisure Centers



Physiotherapy Practises

## Dimensions



### Silverwater

Unit 8, 110-120 Silverwater Rd  
Silverwater NSW 2128  
Ph: 9649 2111

**IMPORTANT** : Always be sure to read your owner's manual and installation guide thoroughly before operation.

[activemobility.com.au](http://activemobility.com.au)

  
Active Mobility

  
ARCHABILITY