



Montar Fold-up Change Table



Montar Fold-Up Change Table is a wall mounted examination table perfect for small spaces. When folded up, it needs very little space and leaves your examination room free for other uses.



Colours



Charcoal

Sizes

Standard

1880mm

Junior

1400mm

Features



Battery Back Up



Weight Capacity

200kg



Safety Rail

Includes a drop down safety rail



Height Adjustable

Adjustable height range of 600mm

Full Features

- The Montar Fold-Up Change Table is a wall mounted examination table perfect for small spaces.
- Needs very little space
- Electrically height adjustable to suit any user
- Upholstery is industrial grade vinyl, which makes it easy to keep clean
- **Complies to IEC60601-2-52**

Optional Features

- Custom Width
- Custom Length
- Custom Vinyl colours
- Wet room upgrade

Specifications

Specifications	Standard	Junior
SWL	200 kg	200kg
Open Overall Length	1880 mm	1400 mm
Open Overall Width	1030 mm	1030 mm
Closed Overall Length	1880 mm	1400 mm
Closed Overall Width	350 mm	350 mm
Bench Surface	1880 x 700 mm	1400 x 700 mm
Height	350-950 mm	350-950 mm
Material	Industrial grade vinyl	Industrial grade vinyl
Colour	Charcoal	Charcoal

Examples



Silverwater

Unit 8, 110-120 Silverwater Rd
Silverwater NSW 2128
Ph: 9649 2111

activemobility.com.au

IMPORTANT : Always be sure to read your owner's manual and installation guide thoroughly before operation.


Active Mobility


ARCHABILITY