



Active Mobility

For more information visit:
activemobility.com.au

Warranty Information
Frame: 10 Years
Molded Foam: 10 Years
Mechanisms: 5 Years
Power: 2 Years
Leather and Fabric: 2 Years

oslö

LUXURY LIFT CHAIRS



SPECIFICATIONS	OSLO ADMIRAL Standard
Maximum Weight Capacity	120kg (SWL)
Overall Width	78cm
Overall Depth	92cm
Overall Height	106cm - 145cm
Seat Width	50cm
Seat Depth	55cm
Seat Height	49cm
Reclined Length	175cm
Unit Weight	61.5kg
SPECIFICATIONS	OSLO ADMIRAL Large
Maximum Weight Capacity	120kg (SWL)
Overall Width	83cm
Overall Depth	92cm
Overall Height	107cm - 147cm
Seat Width	55cm
Seat Depth	55cm
Seat Height	49cm
Reclined Length	175cm
Unit Weight	61.5kg



Oslo Admiral Lift Chair

The Admiral Powered Lift Chair combines refined Scandinavian design and top of the line leather luxury with exclusive patent technology to provide unsurpassed comfort and reliability in the safest and most easily operated lift chair on the market.

A wide range of positioning options can be achieved by the simple hand-held control – choose from up, down and return with capability of adjusting the back and footrest independently. The control also contains a usb port for charging.



Unparalleled Scandinavian design offers superior comfort and function, fitting seamlessly into the décor of any home.



Cold-Cure Moulded Foam encases the interior steel frame and sinuous springs, allowing for body-friendly curvatures, guaranteeing absolute comfort and durability

Exclusive head and neck adjustment system allows your head and neck to rest in an infinite number of comfort positions

A battery backup transformer allows the user to exit the chair even when there is a power outage. Lithium battery pack is also available

Sauvage leather boasts an exceptionally soft feel and elegant looks. Sauvage is premium leather.

Side pocket storage keeps the hand control always within reach and offers ample room for books and magazines

Colours Available*



*Colour may vary in different lighting

With Refined modern looks that are sleek and streamlined, the Admiral is the most luxurious lift chair in the Oslo luxury lift chairs range. Wrapped in Sauvage leather, a naturally milled top grain leather, the Admiral boasts an exceptionally soft feel whilst the unique leather grain provides depth, life, and elegance, guarantying a unique deep and rich look that improves with the age.

With clean lines and meticulous stitching detail, the Admiral Luxury Lift Chair will bring elegance to any home. A timeless design, crafted for ultimate comfort and a refined look, will elevate homes of all styles, from modern to traditional.

Clever Scandinavian design provides for both form and function with limitless options for adjustability from the neck support to the footrest - so you can tailor the chair to your ideal comfort settings. At the touch of a button, you can effortlessly recline, sit and stand for endless comfort, tailored to you.

At the heart of the Admiral and the range of Oslo Luxury Lift Chairs, is safety and ease of operation. Due to a revolutionary patent design that keeps the chair base stable as the seat and back are raised in operation, the Admiral is the safest and most reliable lift chair in the market.

