



## Active Mobility

Unit 8, 110-120 Silverwater Rd  
Silverwater NSW 2128  
Ph: (02) 9649 2111

For more information,  
visit: [activemobility.com.au](http://activemobility.com.au)



SPECIFICATIONS	OSLO ENSIGN Petite
Maximum Weight Capacity	120kg (SWL)
Overall Width	80cm
Overall Depth	87cm
Overall Height	105.5cm - 142.5cm
Seat Width	49cm
Seat Depth	53cm
Seat Height	48cm
Reclined Length	164cm
Unit Weight	64.16kg

SPECIFICATIONS	OSLO ENSIGN Standard
Maximum Weight Capacity	120kg (SWL)
Overall Width	80cm
Overall Depth	90cm
Overall Height	109cm - 144cm
Seat Width	49cm
Seat Depth	53cm
Seat Height	48cm
Reclined Length	164cm
Unit Weight	64.16kg

SPECIFICATIONS	OSLO ENSIGN Large
Maximum Weight Capacity	120kg (SWL)
Overall Width	85cm
Overall Depth	91cm
Overall Height	109 - 144cm
Seat Width	54cm
Seat Depth	53cm
Seat Height	48cm
Reclined Length	164cm
Unit Weight	64.16kg

# oslö

LUXURY LIFT CHAIRS



## Oslö Ensign Lift Chair

The Ensign Powered Lift Chair with semi-lateral back support combines luxurious Scandinavian design and exclusive patent technology to provide unsurpassed comfort and reliability in the safest and most easily operated lift chair on the market.

Unparalleled Scandinavian design offers superior comfort and function, fitting seamlessly into the décor of any home.



Exclusive head and neck adjustment system allows your head and neck to rest in an infinite number of comfort positions

Uniquely designed, contoured back rest with Lumbar and Lateral functions to enhance postural stability

Side pocket storage keeps the hand control always within reach and offers ample room for books and magazines

A battery backup transformer allows the user to exit the chair even when there is a power outage. Lithium battery pack is also available

Colours Available\*



\*Colour may vary in different lighting

A wide range of positioning options can be achieved by the simple hand-held control – choose from up, down and return with capability of adjusting the back and footrest independently. The control also contains a usb port for charging.




Cold-Cure Molded Foam encases the interior steel frame and sinuous springs, allowing significant full body immersion for increased pressure re-distribution whilst guaranteeing absolute comfort and durability

Sleek and streamlined, the clean lines and slender, curved arms of the Ensign Luxury Lift Chair will bring elegance to any home. A timeless design, meticulously finished for a refined look and feel, will suit homes of all styles, from modern to traditional.


Clever Scandinavian design featuring semi-lateral back support, power lumbar and power head rest support allows you to create multi-dimensional sitting and positioning opportunities from head to toe so you can tailor the chair for maximum support and to your ideal comfort settings. At the touch of a button, you can effortlessly recline, sit, and stand for endless comfort, tailored to you.

At the heart of the Ensign and the range of Oslo Luxury Lift Chairs, is safety and ease of operation. Due to a revolutionary patent design that keeps the chair base stable as the seat and back are raised in operation, the Ensign is the safest and most reliable lift chair in the market.


Height Adjustable Head, Neck and Lumbar Support



Lift Chair



Effortless Reclining



Integrated Footstool

