



Warranty Information  
Frame: 10 Years  
Molded Foam: 10 Years  
Mechanisms: 5 Years  
Power: 2 Years  
Leather and Fabric: 2 Years



**Active Mobility**

For more Information visit:  
[activemobility.com.au](http://activemobility.com.au)

SPECIFICATIONS	OSLO MARQUIS
Maximum Weight Capacity	120 kg
Overall Width	79cm
Overall Depth	94 - 98cm
Overall Height	104cm - 154cm
Seat Width	56cm
Seat Depth	57cm
Seat Height	48cm
Reclined Length	176 - 189cm
Wall Clearance (requirement)	43 -51cm
Unit Weight	65.5 kg

# oslö

LUXURY LIFT CHAIRS



## Oslo Marquis Lift Chair

The Marquis Powered Lift Chair with vertical cushioning back support combines luxurious Scandinavian design and exclusive patent technology to provide unsurpassed comfort and reliability in the safest and most easily operated lift chair on the market.

Unparalleled Scandinavian design offers superior comfort and function, fitting seamlessly into the décor of any home.



Manual height and angle adjustable headrest to accommodate a wide array of head support positioning whilst minimizing the risk of inadvertent positioning changes

The backrest features a distinctive, vertically oriented cushioning system that conforms to the natural shape of the spine. This feature ensures optimal alignment and enhanced trunk support, facilitating even distribution of pressure to alleviate discomfort and improve circulation.

Side pocket storage keeps the hand control always within reach and offers ample room for books and magazines

A battery backup transformer allows the user to exit the chair even when there is a power outage. Lithium battery pack is also available

Colours Available\*



Nordic Grey

\*Colour may vary in different lighting

Multi-dimensional positioning with effortless recline, leg elevation and chair lift functions to assist with the user's transfers and repositioning needs.

Easy to use hand hold controller with minimal yet large-sized buttons to increase independence and achieve maximum positioning options for the user



Sleek and streamlined, the clean lines and slender, curved arms of the Marquis Lift Chair will bring elegance to any home. A timeless design, meticulously finished for a refined look and feel, will suit homes of all modern and traditional styles.

Clever Scandinavian design featuring vertical cushioning, manual height and angle adjustable headrest accommodates a wide array of users in one chair and allows you to create multi-dimensional sitting and positioning opportunities from head to toe so you can tailor the chair for maximum support and to your ideal comfort settings. At the touch of a button, you can effortlessly recline, sit, and stand for endless comfort tailored to you.

Safety and ease of operation are at the heart of the Marquis Luxury Lift Chair's design. Its revolutionary patent design keeps the chair base stable as the seat and back are raised in operation, making it the safest and most reliable lift chair on the market.

