



## Active Mobility

Unit 8, 110-120 Silverwater Rd  
Silverwater NSW 2128  
Ph: (02) 9649 2111

For more information,  
visit: [activemobility.com.au](http://activemobility.com.au)

Waterfall  
backrest cushion  
option enhances  
pressure relief  
and individual  
adaptability.



The selection of  
premium materials and  
upholstery adds to its  
refined appearance,  
with choices including  
Mocha and Tuxedo  
engineered fabrics for  
a durable, stylish finish,  
and Anthracite leather  
for a smooth, luxurious,  
high-end aesthetic.

Available in three  
sizes, Petite,  
Standard, and  
Large. This chair  
caters to different  
body proportions,  
ensuring users  
experience precise  
ergonomic fit and  
customized postural  
support.

SPECIFICATIONS	OSLO MINISTER PETITE
Maximum Weight Capacity	150kg
Overall Width	77cm
Overall Depth	92cm
Overall Height	103cm
Backrest Height	72cm
Armrest Height	62cm
Seat Height	46cm
Seat Width	46cm
Seat Depth	46cm
Unit Weight	52 kg
Colours Available	Mocha, Tuxedo, Anthracite

SPECIFICATIONS	OSLO MINISTER STD
Maximum Weight Capacity	150 kg
Overall Width	83cm
Overall Depth	95cm
Overall Height	105cm
Backrest Height	75cm
Armrest Height	65cm
Seat Height	48cm
Seat Width	51cm
Seat Depth	51cm
Unit Weight	55 kg
Colours Available	Mocha, Tuxedo, Anthracite

SPECIFICATIONS	OSLO MINISTER LARGE
Maximum Weight Capacity	180 kg
Overall Width	90cm
Overall Depth	98cm
Overall Height	105cm
Backrest Height	76cm
Armrest Height	68cm
Seat Height	51cm
Seat Width	57cm
Seat Depth	56cm
Unit Weight	58 kg
Colours Available	Mocha, Anthracite

# oslö



## Oslö Minister Adaptable Lift Chair

The Minister  
Powered Lift Chair  
with a choice of  
two comfortable  
cushioning back  
supports combine  
timeless design and  
exclusive technology  
to provide  
unsurpassed  
comfort and  
reliability.

The Oslo Minister Lift Chair is a highly specialised ergonomic seating solution, meticulously developed by a multidisciplinary team including Occupational Therapists, product specialists, and engineers.



Precise  
adjustments can  
be made to the  
headrest

Adaptable  
lateral support  
backrest

Hooks for the  
handset are  
located on  
both sides of  
the chair.

Large storage  
pockets  
to hold  
magazines,  
laptops, books  
etc.

Multi-dimensional positioning with effortless  
recline, leg elevation and chair lift functions  
to assist with the user's transfers and  
repositioning needs.

Colours Available\*



\*Colour may vary in  
different lighting

The Oslo Minister Lift Chair is crafted for  
those who appreciate a compact yet deeply  
ergonomic design.

Engineered for superior posture stability, its  
lateral support backrest adapts effortlessly,  
allowing you to fine-tune both cushioning  
firmness and contour shape.

Powered by four independent motors, this chair  
offers precise adjustments to the headrest,  
lumbar, backrest, and footrest for a fully  
personalised seating experience.

The advanced two-stage lift mechanism ensures  
a seamless, more vertical lift, while the lay-flat  
feature and extended footrest invite ultimate  
relaxation.

Designed for healthy sitting, the built-in lumbar  
support makes every moment in this chair  
effortlessly comfortable.

