

# YOUR OUTDOOR ADVENTURE AWAITS



User-friendly console

Rugged and built to last, the **Outback™** is everything you want in an outdoor scooter! Equipped with front and rear suspension, this aggressive scooter is designed for cruising on walking trails or driving around your property. Enjoy total visibility with complete LED lighting. With a user-friendly console and dual hydraulic brakes for added safety, the Outback is the perfect way to embrace life in the great outdoors!



Active Mobility

(Tel) 02 9649 2111 • Visit [activemobility.com.au](http://activemobility.com.au) for more information.

## 4-Wheel

Model Name	Outback
Weight Capacity	159 kg (350 lbs.)
Maximum Speed <sup>1</sup>	Up to 10 km/h at 159 kg (6.2 mph at 350 lbs.)
Ground Clearance <sup>3</sup>	14 cm (5.5") at motor
Turning Radius <sup>3</sup>	243 cm (94.75")
Overall Length <sup>3</sup>	156 cm (61.25")
Overall Width <sup>3</sup>	83 cm (32.5")
Seat-to-Ground Height	76 - 80 cm (30" -31.5")
Seat-to-Deck Height	50 - 53 cm (19.5"-21")
Front Tyres	37 cm (14.5") tubeless pneumatic
Rear Tyres	37 cm (14.5") tubeless pneumatic
Suspension	Front and rear (CTS Suspension)
Range Per Charge <sup>1,2</sup>	Up to 30 km at 159 kg (Up to 18.6 miles at 350 lbs.) Up to 36 km at 91 kg (Up to 22.6 miles at 200 lbs.)
Standard Seating	Type: Pillow-top, limited recline Weight: 22 kg (48.6 lbs.) Material: Black vinyl Dimensions: Width: 46 cm (18"), Depth: 42 cm (16.5")
Drive System	24-volt DC motor, rear-wheel drive, dual inline
Dual Braking System	Regenerative and electromechanical
Battery Requirements <sup>4</sup>	Size: (2) 100 AH
Battery Charger	Off-board, 8 amp
Warranty	Frame: 2-year limited Drivetrain: 2-year limited Electronics: 1-year limited Battery: 12-months

<sup>1</sup> Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition, and tyre condition. This specification can be subject to a variance of (+ or -) 10%.

<sup>2</sup> Tested in accordance with ANSI/RESNA, MC Vol2, Section 4 & ISO 7176-4 standards. Results derived from theoretical calculations based on battery specifications and drive system performance. Test conducted at maximum weight capacity.

<sup>3</sup> Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of (+ or -) 3%.

<sup>4</sup> AGM or gel cell type required.

NOTE: All specifications subject to change without notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tyre pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3 %. Special order seats or non-standard dimensions will not be fold-flat and depth-adjustable.

Theoretical range is calculated under test conditions in accordance with European standards. Actual range in normal use will depend on many factors, including the condition of the vehicle and its batteries, the weight of the driver, correct tyre pressures, ambient temperature, and the gradient and surface of the road or pavement. Dimensions listed are for power base only. Overall measurements will vary based on seating and accessory selections.

## Key Features

- Limited recline, black vinyl, high-back seat with lumbar and pillow-top headrest, flip-up height, and angle-adjustable armrests
- CTS Suspension includes adjustable shocks for greater comfort
- Backup sensors with LED light
- Dual hydraulic brakes

## Additional Features

- User-friendly LED console displays time, temperature, kilometres driven, and trips taken
- Adjustable delta tiller with ergonomic, wraparound handles / storage compartment included
- USB charger built into the tiller for convenient charging of smart phones and portable devices
- Durable and stylish rear bumper
- Easily accessible tie-down points (for transport of unoccupied scooter)
- Full LED lighting package includes headlights, hazard lights, rear backup sensor LED lights, and directional signals
- XLR charger port

## Available Colours



Stone  
Grey

## Accessories

- Rear bag / walking stick holder
- Oxygen tank holder\*
- Mobile phone holder
- Flag
- Walker holder\*
- Cup holder
- Saddlebag

\* Must have Pride rear bag accessory to attach



Active Mobility Systems Pty Ltd  
Unit 8/110-120 Silverwater Road  
Silverwater NSW 2128