



PRIDE® Pathrider 140 XL

An exciting blend of power, precision and exceptional style, the Pathrider 140 XL is the perfect outdoor scooter.




Colours Available



Yellow



Candy
Apple Red



Silver

Features



Seat
High back
captain seat
with pillow top



Wheels
4 Wheel



Weight Capacity
181kg
Bariatric suitable



Lighting
LED headlights,
front shield, indicators,
hazard lights

Full Features

- Powerful drivetrain with 24-Volt DC motor
- Full suspension for optimum handling
- Sealed hydraulic brakes
- Large 33cm pneumatic tyres for excellent outdoor performance
- Easy drive tiller with wraparound handles
- Easily accessible tie down points (for transport of unoccupied scooter)
- Rearview mirror
- Front shield
- High back pillow top seat with headrest
- Front & rear suspension
- Front basket
- LED lower headlight
- Hazard lights / indicators
- Sleek, sporty style

Specifications

SPECIFICATIONS	PRIDE PATHRIDER 140XL
Maximum Weight Capacity	181kg
Maximum Speed ¹	15km/h
Maximum Range ^{1,4}	Up to 35km with 75AH batteries
Length ²	142.8cm
Width ²	69.2cm
Scooter Weight (without battery) ^{2,5}	117.70kg
Ground Clearance ²	13.33cm
Seat to Ground	72.39cm - 68.58cm
Turning Radius ²	Right: 210.8cm Left: 208.3cm
Battery Requirements ³	2 x 12-Volt, Deep Cycle 75AH Batteries
Charger	8Amp, Off Board
Seating	Pillow Top, 21.77kg, Black Vinyl, Width 49.5cm, Depth 44.5cm
Motor	Rear Wheel Drive, Sealed Transaxle, 24-volt DC Motor

¹Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition and tyre condition. These specifications can be subject to a variance of (+/- 10%).

² Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of (+ or - 3%).

³AGM or gel cell type required.

⁴Tested in accordance with ANSI/RESNA, WC Vol2, section 4 & ISO 7176-4 standards. Results derived from theoretical calculations based on battery specifications and drive system performance.

Test conducted at maximum weight capacity.

⁵Battery weight may vary based on manufacturer.

IMPORTANT : Always be sure to read your owner's manual thoroughly before operating your machine.

WARNING : Do not exceed 10km/ph on public roads or paths.

Silverwater

Unit 8, 110-120 Silverwater Rd
Silverwater NSW 2128
Ph: 9649 2111

Newcastle

7 Poynton Pl
Thornton, NSW 2322
Ph: 4047 3608

activemobility.com.au



PRIDE[®]