



Active Mobility

SKU 1559017

Quest

Ready for transport and easily folded, the 4-wheel Quest Folding Scooter is perfect for the active lifestyle.



Colours



White

Features



Speed
6 km/h



Weight Capacity
115kg



LED Lighting
Front and rear LED lighting



Portable
This scooter is able to be folded for easy transportation

Full Features

- great performance indoors and out at a maximum speed of 6 km/h
- Safe working load of 115kg
- Range per charge of up to 20 km
- A great space saver that adjusts to fit perfectly inside any wardrobe or vehicle boot
- Front and rear LED lighting
- Unique folding design that allows the scooter to be folded in three easy steps
- Adjustable tiller with easy to grip adjustment knob
- Can be stowed upright or layed down
- Articulated front end adds greater stability

Specifications

SPECIFICATIONS	iRIDE
Weight Capacity	115kg
Maximum Speed	6km/h
Ground Clearance	6 cm
Turning Radius	82cm
Overall Length	93cm
Overall Width	48.5cm
Front Tyres	4 x 18 cm (1.6" x 7") solid
Rear Tyres	5 x 20 cm (2" x 8") solid
Range per Charge	Up to 20km
Total weight without batteries	23.5kg
Battery Requirements	(1)24-volt (Lithium-ion) Size: (1) 11.5 AH (Lithium-ion) Weight: 11.5 AH (Lithium-ion) Each 2.3 kg (5 lbs.) (One Required)
Battery Charger	36V lithium, off-board

Warranty

Warranty Information	
Scooter	12 Month Warranty for all electronic parts incl controllers & battery chargers
Scooter Frames	24 Month warranty on the frame incl platform , fork , seat post , drivetrain, differential ,motor , brakes and welds.
Batteries	6 Month Warranty

IMPORTANT : Always be sure to read your owner's manual thoroughly before operation.

Silverwater
 Unit 8, 110-120 Silverwater Rd
 Silverwater NSW 2128
 Ph: 9649 2111

Newcastle
 7 Poynton Pl
 Thornton, NSW 2322
 Ph: 4047 3608

activemobility.com.au

