



Hit the streets in style with the Pride Sportrider. The Sportrider is a unique 3-wheel motorcycle inspired, outdoor, sporty scooter




## Colours



Black

## Features



**Suspension**  
Front and rear suspension



**Weight Capacity**  
180kg



**Front Lighting**  
Powerful front and rear lighting kit



**Storage**  
Comes included with a large rear storage pod

# Full Features

- A full digital display that includes a speedometer/odometer
- Safe working load of 180kg
- Equipped with front and rear suspension, the Sportsrider can comfortably cover most terrain
- Enjoy total visibility with a powerful front and rear lighting kit
- High-back, adjustable seat
- Large rear storage pod

# Specifications

SPECIFICATIONS	OUTBACK
Weight Capacity	180kg
Maximum Speed	10km/h
Ground Clearance	13 cm (5.1")
Turning Radius	157.3cm
Overall Length	165cm
Overall Width	78cm
Seat to Floor Height	55cm
Front Tyres	35.6cm (14.") pneumatic
Rear Tyres	35.6cm (14.") pneumatic
Suspension	Front and Rear
Range per Charge	Up to 45-50km
Total Unit weight with Batteries	188.5 kg
Seating	Weight: 22.5kg /Material: Black Vinyl/Dimensions: 50.4 x 50 cm
Motor Output	1300W 4-pole 3.5 HP
Brakes	Regenerative and electromechanical
Battery Requirements	MK Gel Batteries, Size: 70 AH, Weight 34kg each
Battery Charger	Off-board, 8 amp
Incline Capabilities	10°

# Warranty

Warranty Information	
Scooter	12 Month Warranty for all electronic parts incl controllers & battery chargers
Scooter Frame	24 Month warranty on the frame incl platform , fork , seat post , drivetrain, differential ,motor , brakes and welds.
Battery	6 months

**IMPORTANT :** Always be sure to read your owner's manual thoroughly before operation.

**Silverwater**  
Unit 8, 110-120 Silverwater Rd  
Silverwater NSW 2128  
Ph: 9649 2111

**Newcastle**  
7 Poynton Pl  
Thornton, NSW 2322  
Ph: 4047 3608

[activemobility.com.au](http://activemobility.com.au)

