




# iRIDE

Journey far and wide with the iRIDE. This lightweight scooter is super portable and perfect for travelling.



## Colours



Arctic Ice



Raspberry

## Features



**Light Weight**  
14.8 kg<sup>1</sup>



**Weight Capacity**  
113kg



**Hand Brake**  
Equipped with a manual hand brake



**Portable**  
This scooter is able to be folded for easy transportation

# Full Features

- At only 14.8 kg (without seat/batteries), this lightweight scooter is super portable and perfect for travelling.
- Safe working load of 113kg
- Equipped with both a manual hand brake and supplemental foot brake, the iRIDE delivers a safe and comfortable ride
- Digital LCD display
- Foldable tiller with easy to grip adjustment knob
- Airline-safe lithium-ion battery pack
- Height adjustable seat
- Up to 6.2km/h

# Specifications

SPECIFICATIONS	iRIDE
Weight Capacity	113kg
Unit Weight <sup>1</sup>	14.8kg
Maximum Speed	6.2km/h
Ground Clearance	3.2 cm (1.25") at rear disc
Turning Radius	85cm
Overall Length	96.5cm
Overall Width	48.8cm
Seat to Ground Height	45.7-50.8cm
Seat to Deck Height	25.2-30.2cm
Front Tyres	16.5cm (6.5") solid
Rear Tyres	17.1cm (6.75") solid
Suspension	N/A
Range per Charge	Up to 13.3km (at 113.3kg)
Total weight without batteries	19.38kg
Weight of Heaviest Peice	14.8kg (base only)
Battery Requirements	Size: (1) 36V 5.2 AH lithium ion Weight: 2.14 kg (4.7 lbs.)
Dual Braking System	Manual, dual disc
Battery Charger	36V lithium, off-board

# Warranty

Warranty Information	
Frame	2 year limited
Drivetrain	1 year limited
Electronics	1 year limited
Battery	6 months

<sup>1</sup> Weight is calculated with seat and batteries removed.

**IMPORTANT :** Always be sure to read your owner's manual thoroughly before operation.

**Silverwater**  
Unit 8, 110-120 Silverwater Rd  
Silverwater NSW 2128  
Ph: 9649 2111

**Newcastle**  
7 Poynton Pl  
Thornton, NSW 2322  
Ph: 4047 3608