




rolio^{DUO}

Ergonomically designed with sturdy lightweight construction and aluminum frame, the premium Rolio Duo 2-in-1 Walker is able to convert into a wheelchair in seconds and serve multiple purposes in one unit.



Colours Available



Rose Gold



Olive Grey

Features



2-In-1
The Rolio Duo is able to convert from a walker to a wheelchair for a seamless transition from standing to sitting



Weight Capacity
136kg



Walker Weight
9.9kg

Full Features

- Premium 2-in-1 walker and wheelchair ergonomically designed to serve multiple purposes in one unit
- One touch, easy-fold and unfold mechanism for convenient transport and storage
- Safety breaks and push down locks for additional safety and stability when rising from a seated position
- Sturdy lightweight construction, aluminum frame weighing under 10kg
- Lift up swing away locking legrests with lightweight flip-up footplates
- Patent swivel backrest allowing maximum back support in both walker & transit chair positions
- Height adjustable handles allowing up to 130mm of additional handle height
- Large removable tote bag with zip allows for ample storage of a handbag and belongings
- Kerb handler feature enables easy manouvering over kerbs, steps and uneven ground

Specifications

SPECIFICATIONS	AZURE ROLIO DUO
Maximum Weight Capacity	136kg
Width with Legrest Folded	56cm
Width with Legrest Extended	61cm
Depth with Legrest Folded	74cm
Depth with Legrest Extended	99cm
Height	85.3mm
Unit Weight	9.9kg
Folded Width with Legrest Folded	28cm
Folded Width with Legrest Extended	29cm
Folded Depth with Legrest Folded	74cm
Folded Depth with Legrest Extended	99cm
Folded Height	85.3cm
Handle Height	84.5cm - 97cm
Transfer Handle Height	79cm
Seat Width	46cm
Seat Depth	22.5cm
Seat Height	53cm
Front Wheel	25cm / 10"
Rear Wheel	20cm / 8"

IMPORTANT : Always be sure to read your owner's manual thoroughly before operation

WARRANTY INFORMATION: Rolio is issued with a 12 month product warranty

Unit 9, 110-120 Silverwater Rd
Silverwater NSW 2128
Ph: 9649 2111

azuredesigns.com.au

Azure | **rolio**^{DUO}
DESIGNS