



rollo

Unit 8, 110-120 Silverwater Rd
Silverwater NSW 2128
Ph: (02) 9649 2111

For more information,
visit: activemobility.com.au

rollo TALL

15878640 - Graphite Grey



Active Mobility



rollo TALL

Designed for function and style

Boasting ergonomic design, supreme adjustability, comfort, and class, rollo tall is the perfect combination of function and style in mobility.

Rolio Tall's **ergonomic shape and taller handles** fit perfectly in the hand delivering a softer and smoother handling. Glide over most terrains with total confidence and ease!

Lightweight aluminium frame allows rollo to be easily transported anywhere making it a trusty companion for shopping trips, walks in the parks or on holiday.

Ultra-comfort seat crafted from high quality durable materials provides supreme cushioning comfort during rest breaks.



Rolio's **handy tote bag** provides ample space for a handbag and essentials so you can travel light and carefree wherever you go.



Designed for function and style

The **rollo walker makes a statement** with a sleek, streamlined frame available in the new Tall model, to accommodate for taller users. The Rollo tall has a seat to floor height of 60cm (compared to Rollo standard : 54cm), and features a higher handle height of 93cm which can be extended to 105cm.



Patent adjustable padded backrest providing up to 80mm of seat depth adjustment and 160mm height adjustment allowing the user to customise rollo to their ideal comfort settings.



Simple one touch folding and unfolding mechanism takes away any fuss during travel, allowing you to breeze through your day with ease.



136kg
Max User Weight