




Active Mobility

SKU 15878610




rollo

Boasting ergonomic design, supreme adjustability, comfort and class, the rollo walker is the perfect combination of function and style in mobility.



Colours Available



Rose Gold



Olive Grey

Features



Adjustable Backrest

Up to 80mm of seat depth and 160mm height adjustment



Weight Capacity

136kg



Walker Weight

7.67kg

Full Features

- Sturdy construction lightweight, aluminum frame weighing under 8kg
- One touch, easy-fold and unfold mechanism for convenient transport and storage
- Safety breaks and push down locks for additional safety and stability when rising from a seated position
- Ergonomically designed hand grips with safety reflectors
- Ultra-comfort seat with supreme cushioning
- Patent adjustable padded backrest providing up to 80mm of seat depth adjustment and 160mm height adjustment
- Height adjustable handles allowing up to 130mm of additional handle height
- Large removable tote bag with zip allows for ample storage of a handbag and belongings
- Kerb handler feature enables easy manouvering over kerbs, steps and uneven ground

Specifications

SPECIFICATIONS	AZURE ROLIO
Maximum Weight Capacity	136kg
Width	56cm
Depth	74cm
Height	85.3mm
Unit Weight	7.67kg
Folded Width	28cm
Folded Depth	74cm
Folded Height	85.3cm
Handle Height	84.5cm - 97cm
Seat Width	46cm
Seat Depth	22.5cm
Seat Height	53cm
Front Wheel	25cm / 10"
Rear Wheel	20cm / 8"

IMPORTANT : Always be sure to read your owner's manual thoroughly before operation

Warranty: 12 Month Warranty

Unit 8, 110-120 Silverwater Rd
Silverwater NSW 2128
Ph: 02 9649 2111

7 Poynton Place,
Thornton NSW 2322

activemobility.com.au


Active Mobility

