



Active Mobility

SKU 15878630



# rollo+

Boasting ergonomic design, supreme adjustability, comfort and class, the rollo+ walker is the perfect combination of function and style in mobility.



## Colours Available



Olive Grey

## Features



### Adjustable Backrest

Up to 80mm of seat depth and 160mm height adjustment



### Weight Capacity

180kg



### Walker Weight

9.45kg

# Full Features

- Suitable for
- Sturdy construction lightweight, aluminum frame weighing under 10kg
- One touch, easy-fold and unfold mechanism for convenient transport and storage
- Safety breaks and push down locks for additional safety and stability when rising from a seated position
- Ergonomically designed hand grips with safety reflectors
- Ultra-comfort seat with supreme cushioning
- Patent adjustable padded backrest providing up to 80mm of seat depth adjustment and 160mm height adjustment
- Height adjustable handles allowing up to 130mm of additional handle height
- Large removable tote bag with zip allows for ample storage of a handbag and belongings
- Easy manouverability and smooth handling over most terrains

# Specifications

SPECIFICATIONS	AZURE ROLIO+
Maximum Weight Capacity	180kg
Width	66.5cm
Depth	73cm
Height	92cm
Unit Weight	9.45kg
Folded Width	28cm
Folded Depth	73cm
Folded Height	92cm
Handle Height	92cm - 105cm
Seat Width	56cm
Seat Depth	30.5cm
Seat Height	60cm
Front Wheel	25cm / 10"
Rear Wheel	20cm / 8"

**IMPORTANT :** Always be sure to read your owner's manual thoroughly before operation

**Warranty:** 12 Month Warranty

Unit 8, 110-120 Silverwater Rd  
Silverwater NSW 2128  
Ph: 9649 2111

[activemobility.com.au](http://activemobility.com.au)

  
Active Mobility

