



NEW

rolio CARBON ULTIMATE


With its ergonomic, lightweight design, the **Rolio Carbon Ultimate Walker** combines advanced carbon fibre technology with luxurious leather detailing, allowing you to move with confidence, fold with ease, and rest comfortably no matter where the day takes you.



Colour



White Frame



Brown Leather

Features



Weight Capacity
150kg



Walker Weight
6.5kg



Removable Bag
An easy slide on and off handle that securely holds your bag



Large Castors
Castors for smoother steering and obstacle clearance

The **Rolio Carbon Ultimate** is the pinnacle of walking aids, pairing an ultra-light carbon fibre frame with exquisitely finished leather detailing for a truly premium experience. Designed for those who value refinement as much as performance, it moves effortlessly, folds seamlessly and offers a beautifully comfortable place to rest. With smooth handling over uneven ground and thoughtful details like a discreet zip bag, the Rolio Carbon Ultimate delivers confidence, comfort and unmistakable sophistication, supporting your independence with style, and backed by a 12-month warranty.

Full Features

- Strong, lightweight carbon fibre frame weighing only 6.5kg
- One touch, easy-fold and unfold mechanism for convenient transport and storage
- Safety breaks and push down locks for additional safety and stability when rising from a seated position
- Ergonomically designed hand grips with safety reflectors
- Ultra-comfort leather seating with supreme cushioning
- Patent adjustable padded backrest providing up to 160mm height adjustment
- Height adjustable handles allowing up to 130mm of additional handle height
- Large removable leather bag with zip allows for ample storage of a handbag and belongings
- Kerb handler feature enables easy manoeuvring over kerbs, steps and uneven ground

Specifications

SPECIFICATIONS	AZURE ROLIO CARBON ULTIMATE
Maximum Weight Capacity	150kg
Width	56cm
Depth	74cm
Unit Weight	6.5kg
Folded Width	28cm
Folded Depth	74cm
Overall Height	84.5cm - 97cm
Seat Width	46cm
Seat Depth	22.5cm
Seat Height	53cm
Front Wheel	25cm / 10"
Rear Wheel	20cm / 8"
Material	Carbon Fibre
Colour	White

IMPORTANT : Always be sure to read your owner's manual thoroughly before operation

Unit 9, 110-120 Silverwater Road
Silverwater NSW 2128

Ph: 02 9737 0144

azuredesigns.com.au

sales@azuredesigns.com.au

Azure
DESIGNS

rolio CARBON
ULTIMATE