

PORTABLE ROLLATOR WALKER

ROVENA

ROVENA offers unmatched **comfort, stability, and convenience.**

Whether you're running daily errands or enjoying a leisurely stroll outdoors, its **durable steel frame** ensures **smooth and effortless** steering.



Colours Available



Blue



Dark Red



Black

Features



Fold Design

The Rovena incorporates a compact front fold design that makes it easy to store and transport



Weight Capacity
150kg



Integrated Storage
Comes with a built-in storage bag conveniently tucked beneath the walker

Full Features

- Sturdy steel frame, weighing under 9kg
- Durable, smooth-rolling 8-inch castors designed for both indoor and outdoor use
- Push-to-lock brakes for added safety — just push down to lock the wheels in place when stationary
- Thicker chair cushion ensures a more comfortable sitting experience, ideal for longer periods of use
- Thicker backrest for comfort & support
- Adjustable handle height from 840mm to 980mm, improving ergonomics and reducing strain
- Available in three seat widths — 41 cm, 46 cm, or 51 cm — to suit a wide range of body types and user preferences
- Integrated baggage underneath the seat with a max load of 5kg (corrected spelling from “Intergrated” to “Integrated”)
- Front fold-open mechanism allows for quick and compact storage
- Integrated fold-open chair provides access to storage and offers a reliable resting spot

SPECIFICATIONS	ROVENA
Material	Steel
Max. load	150kg
Max. load of the bag	5 kg
Total weight without shopping bag	8.7
Shopping	0.6kg
Seat width	450mm
Height of hand grips (adjustable)	840-980mm
Height of seat	530mm
Brakes	2 x maneuver brakes with parking function
Wheels (diamter x width)	front: 200 x 36 mm rear: 200 x 36 mm
Unfold size	660*620*840 mm
Fold size	910*620*220 mm
Package size	610*203*825 mm
Loading QTY	609 PCS/HC
Total weight with carton box	10.8kg/box
Standard configuration	One shopping bag and one seat

Unit 9
110-120 Silverwater Rd
Silverater NSW 2128
02 9737 0144

Azure
DESIGNS

azuredesigns.com.au sales@
azuredesigns.com.au

ROVENA

