




## 889SL

The SHOPRIDER™ 889SL is a premium model from the mid size scooters available





### Colours




### Features

 **Speed**  
10 km/h

 **Weight Capacity**  
136kg

 **LED Lighting**  
Front and rear LED lighting

 **Incline Capability**  
Up to 9° (91kg)  
Up to 6° (113kg)

# Full Features

- Great performance outdoors at a maximum speed of 10 km/h
- Safe working load of 136kg
- Range per charge of up to 40 km
- Front and rear LED lighting, including hazard lights
- Large Captain's Seat, that can slide and rotate for easy transfer
- Battery level warning light
- Front and rear bumper pads
- Front and rear suspension
- Adjustable width and foldable arm-rests
- High/low speed switch
- Anti-tip wheels

# Specifications

SPECIFICATIONS	889SL
Weight Capacity	136kg
Maximum Speed	10km/h
Length	1300mm
Height	1200mm
Safe Climbing Angle	9°
Max Climbing Angle	12°
Charger	4 Amp
Motor Size	1.20Hp (4 pole)
Battery Capacity	40Amp (minimum Amp)
Maximum Range	up to 40 km
Wheel Size	330mm
Turning Radius	1440mm
Weight (with battery)	99kg

# Warranty

Warranty Information	
Scooter Frames	10 Years Structural Frame 18 Months Drive Train and Batteries
Batteries (Gel/Hybrid)	18 Months
Scooter Parts and Accessories	1 Year

**IMPORTANT :** Always be sure to read your owner's manual thoroughly before operation.

**Silverwater**  
Unit 8, 110-120 Silverwater Rd  
Silverwater NSW 2128  
Ph: 9649 2111

**Newcastle**  
7 Poynton Pl  
Thornton, NSW 2322  
Ph: 4047 3608

[activemobility.com.au](http://activemobility.com.au)

