




# SHOPRIDER® AllRounder 889


Jam packed with features, the AllRounder has everything you want in a mobility scooter.



## Colours Available



Blue



Burgundy Red



Champagne Silver

## Features



**Seat**  
High back contoured captain seat with headrest



**Weight Capacity**  
159kg  
Bariatric suitable



**Wheels**  
4 Wheel



**Lighting**  
LED headlights, rear lights, indicators, hazard lights


# Full Features

- Powerful 4 pole motor
- High back contoured captain's seat with headrest
- Adjustable steering tiller
- Adjustable width arm rests
- Sliding rotating height adjustable seat
- Fold arm rests
- Variable seat height
- Front shopping basket
- Horn
- LED battery level indicator
- Front chrome bumper
- Front & rear adjustable suspension
- Front & rear lights
- Hazard lights / indicators
- High / low speed switch
- Anti tip wheels

# Specifications

SPECIFICATIONS	SHOPRIDER ALLROUNDER 889
Maximum Weight Capacity <sup>1</sup>	159kg
Maximum Speed	10km/h
Maximum Range <sup>2</sup>	Up to 40km
Length	130cm
Width	64cm
Scooter Weight (with battery)	99kg
Turning Radius	144cm
Battery Requirements	2 x 40Amp
Charger	8Amp
Seating	High Back Air Flow Seat,
Motor Size	1.5Hp (4 pole)

<sup>1</sup>Single User



<sup>2</sup>Driving range varies with the factors shown below:

- (a) The weight of the occupant
- (b) Level ground & ambient temperature
- (c) Battery condition / Type of charger
- (d) Driving Style

### Warranty Information:

Scooters: 10 year structural frame, 18 months drivetrain and batteries

Parts and Accessories: 12 months

Batteries (Gel/Hybrid): 18 months

**IMPORTANT:** Always be sure to read your owner's manual thoroughly before operating your machine.

**WARNING:** Do not exceed 10km/ph on public roads or paths.

### Silverwater

Unit 8, 110-120 Silverwater Rd  
Silverwater NSW 2128  
Ph: 9649 2111

### Newcastle

7 Poynton Pl  
Thornton, NSW 2322  
Ph: 4047 3608

[activemobility.com.au](http://activemobility.com.au)



**SHOPRIDER**