



Active Mobility



BURGUNDY RED

BLUE

ROCKY 4

1569007

FEATURES

- High back contoured seat with headrest
- Gas adjustable steering tiller
- Loop handle with finger operated controls
- High / Low speed settings
- Rear wheel drive
- Electromagnetic brake
- Front + Rear (Adjustable) suspension
- Twin Headlights (LED)
- Rear lights
- Indicators
- Hazard lights
- Large battery capacity
- Battery power indicator
- USB power outlet
- Sliding / Rotating / Height adjustable seat
- Lift up & adjustable width arm rests
- Removable shopping basket

SPECIFICATIONS

Overall: 1430mm (L) x 680mm (W) x 1300mm (H)

Turning Circle: 1800mm

Speed: 10km/h

Maximum range: Up to 40km

Safe climb angle: 9°

Maximum climb angle: 12°

Motor: 2.0Hp high torque motor (4 Pole)

Batteries: 2 x 75 Amp minimum (supplied)

Charger: 8 Amp

Tyres (Pneumatic): Ø330mm - Front + Rear

Product Weight: 146kg - with batteries

Maximum User Weight: 225kg - on flat ground



0° - 225kg
6° - 180kg
9° - 159kg
12° - 136kg